

# BLOOD PRESSURE

## FITTINGS



BPL Series 1/4 Turn Cuff Connectors

**PG-37**



BPF Series Female Bayonet Style Quick Connects

**PG-37**



Bayonet Male Cuff Connectors

**PG-38**



Bayonet Panel Mount Connectors

**PG-38**



Dinamap® Monitor Adapters

**PG-38**



Dinamap® Cuff Connectors

**PG-39**



## **BPL** Series

The BPL Series quick connect fittings offers superior ergonomics - a requirement for many designers in their effort to accommodate end users. The BPL's ease-of-use stems largely from a winged male fitting and Value Plastics' patented 10-point grip, making the quarter-turn connection simple and easy to accomplish. A solid snap, that can be both heard and felt, will then reassure the end user that a secure fit has been made.

Another unique characteristic of the BPL Series is the true tube stop on the rotating half of the two connectors. This allows tubing to be fully assembled onto the barb without its edge pushing against the body of the connector.

### Ergonomic wing-grip design

- Easier to connect and disconnect than market similars
- 360° of tubing rotation

### Quarter-turn design with closure feedback

- Positive reinforcement of secure connection

### Sizes to fit 1/8" (3.2 mm) and 5/32" (4.0 mm) ID tubing

- Fits almost all cuffs
- 5/32" hose barb perfectly fits most common cuff tube size

### 200 Series single barbs

- Ease of assembly
- Enhanced resistance to tube pulloff

### Medical grade ABS material (fluid path)

- Withstands bulk sterilization processes
- Holds up under demanding application conditions

### Acetal (Wing-grip shroud)

- Added strength