

Material - Manufacturer & Grade		VP Code	Regulatory Compliance of Base Resin*		RoHS Compliant	BPA Free	DEHP Free	Animal Derivative Free*
ABS (Lace White) - INEOS Lustran 348-000000	-81		USP Class VI Certified; ISO 10993 Part I Compliant	Yes	Yes	Yes	No	
ABS, MABS (Methyl Methacrylate) - BASF Terlux 2802 TR	-8003, SCV Series		USP Class VI Certified	Yes	Yes	Yes	Yes	
ABS (Sno White) - INEOS Lustran 348-012002	-8012		USP Class VI Certified; ISO 10993 Part I Compliant	Yes	Yes	Yes	No	
Acetal - DuPont Delrin 500P	-10		21 CFR 177.2480 FDA Approved	Yes	Yes	Yes	No	
Acrylic Copolymer - BASF Ultraclear N 23020 003	-1006, -1007		21 CFR 177.2470 FDA Approved	Yes	Yes	Yes	No	
Acrylic - Cyro Cyrolite Med 2	-50		USP Class VI Certified; ISO 10993 Compliant; 21 CFR 177.1010 FDA Approved	Yes	Yes	Yes	Yes	
Acrylic (PA 66) - DuPont Zytel 101F	-0, -1, -2, -3, -4, -5, -06, -07, -08B, -8, -8B, -CM030		21 CFR 177.1500 FDA Approved	Yes	Yes	Yes	No	
Acrylic (PA 66) - DuPont Zytel SC310 NC010	-U1, -U5, -U01, -U05		USP Class VI Certified; ISO 10993 Compliant; 21 CFR 177.1500 FDA Approved	Yes	Yes	Yes	No	
Acrylic (Glass Reinforced) - DuPont Zytel IGFES171	-U71, -U75		USP Class VI Certified; ISO 10993 Compliant; 21 CFR 177.1500 FDA Approved	Yes	Yes	Yes	No	
Acrylic (Glass Reinforced) - DuPont Zytel 70G33L	-7, -72		21 CFR 177.1500 FDA Approved	Yes	Yes	Yes	No	
Acrylonitrile (Clear) - Bayer Makrolon 2558 550115	-9		USP Class VI Certified; ISO 10993 Compliant; 21 CFR 177.1580 FDA Approved	Yes	No	Yes	No	
Acrylonitrile (Clear) - Bayer Makrolon Rx 1805 451118	-9010, SCV Series		USP Class VI Certified; ISO 10993 Compliant	Yes	No	Yes	Yes	
Acrylonitrile (Tinted) - Styron Calibre Megarad 2081-15	-9002, DCV Series		ISO 10993 Compliant	Yes	No	Yes	No	
Acrylonitrile - LyondellBasell Purell PE 1840 H	SCV Series		Please call for details	Yes	Yes	Yes	Yes	
Acrylonitrile - Flint Hills P5M6K-080	-6005, -6006		USP Class VI Certified; 21 CFR 177.1520 FDA Approved	Yes	Yes	Yes	Yes	
Polysulfone - Solvay Udel P1700 CL2611	-40		USP Class VI Certified; ISO 10993 Compliant; 21 CFR 177.1655 FDA Approved	Yes	Yes	Yes	Yes	
PVDF - Arkema Kynar 1000 HD	-J1A		USP Class VI Certified; 21 CFR 177.2510 FDA Approved	Yes	Yes	Yes	Yes	
SAN - BASF 388S TR 77770	SCV Series		ISO 10993 Compliant	Yes	Yes	Yes	No	
Silicone - Wacker Elastosil LR 3003	SCV Series		USP Class VI Certified	Yes	Yes	Yes	Yes	

*As determined by resin manufacturer