







VPS5401068N	NUMBER	MATERIAL
	VPS5401068N	Polycarbonate with Silicone Disk
Check Valve, Male/Female Luer, 420 mL/min flow rate at 1 M head 30 psi back pressure, 1/4 psi inlet cracking pressure		

VPS5401020N	NUMBER	MATERIAL
	VPS5401020N	Polycarbonate with Silicone Disk
Check Valve, Male/Female Luer, 129 mL/min flow rate at 100mmHg, 45 psi back pressure, 1-5 psi inlet cracking pressure		


VPS5401036SN	NUMBER	MATERIAL
	VPS5401036SN	Polycarbonate with Silicone Disk
Reflux Check Valve, Male/Female Luer, 350mL/min flow rate at 1 M head, 45 psi back pressure, opens at 284 mmHg		


VPS5401001N	NUMBER	MATERIAL
	VPS5401001N	ABS Housing with Silicone Disk
Check Valve, Male Slip Luer, 420 mL/min flow rate at 1M head, 30 psi back pressure, 1/4 psi inlet cracking pressure		


VPM2930601N	NUMBER	MATERIAL
	VPM2930601N	White ABS
Tube Clamp with split, 1/8" (3.2 mm) OD Tubing		


VPS3901074N	NUMBER	MATERIAL
	VPS3901074N	Polyisoprene with ABS Acrylic Body
IN-4000 Intermittent Injection Site to Male Luer Lock		


VPM1060200N	NUMBER	MATERIAL
	VPM1060200N	Clear Low Density Polyethylene
Long Spike Guard Closed at Tip		


VPM0480201N	NUMBER	MATERIAL
	VPM0480201N	White Polyethylene
Luer Tip Syringe Cap		

VPB0850400N	NUMBER	MATERIAL
	VPB0850400N	Clear ABS
Female Luer Thread Style Coupler		


VPM0031400N	NUMBER	MATERIAL
	VPM0031400N	Clear Acrylic
Female Luer Lug Style, Bond-in 5/32" ID		


VPM3341750SN	NUMBER	MATERIAL
	VPM3341750SN	Blue Rad. Stable Polycarbonate
Tri-female Luer Y Connector to Male Luer Slip		

VPM1080303N	NUMBER	MATERIAL
	VPM1080303N	Blue Polyethylene
Male Luer Integral Lock Ring Plug		

VPB1000055N	NUMBER	MATERIAL
	VPB1000055N	Polycarbonate Body with Polyethylene Handle
One way Stopcock, Male Slip Luer to Capped Female Luer Thread Style		

VP455980	NUMBER	MATERIAL
	VP455980	Polycarbonate Body with Polyethylene Handle
One Way Stopcock, Capped Male Luer Lock to Capped Female Luer Thread Style *Sterile Version*		

VPB1000079N	NUMBER	MATERIAL
	VPB1000079N	Polycarbonate Body with Polyethylene Handle
Three Way Stopcock, (2) Capped Female Luer Thread Style and Capped Male Luer Lock		

VPB15311LG2	NUMBER	MATERIAL
	VPB15311LG2	Rad. Stable Polycarbonate Body with White HDPE Handle
3-Way Stopcock, (2) Female Luer Thread Style & (1) Male Luer with Rotating Lock Ring (NO caps)		

Australian Distributors: **HRMalytic** www.chromtech.net.au sales@chromtech.net.au
 A&M 74 643 445 058 **PTY LTD** Tel: (03) 9762 2034 Fax: +61 3 9761 1169
TECH nology ... for ALL Your chromatography supplies !