

Material Suffix Key

-001 Black Nylon over White Nylon
-81 Lace White ABS

-1 White Nylon
-8012 White ABS

-5 Blue Nylon

-10 Acetal

BPM013 (-1)

Blood Pressure Monitor Bayonet Connector to 500 Series Barb, 1/8" (3.2 mm) ID Tubing



works with...

**BPM240 (-1)**

Blood Pressure Monitor Bayonet Connector to 200 Series Barb, 5/32" (4.0 mm) ID Tubing



works with...

**BPM025 (-1)**

Blood Pressure Monitor Bayonet Connector to 500 Series Barb, 5/32" (4.0 mm) ID Tubing



works with...

**BPM035 (-1)**

Blood Pressure Monitor Bayonet Connector to 500 Series Barb, 3/16" (4.8 mm) ID Tubing



works with...

**BPMMTLL (-1 -5)**


Blood Pressure Monitor Bayonet Connector to Male Luer Integral Lock Ring



works with...

**BPAM (-001)**


Blood Pressure Monitor Bayonet Connector to Internal 5/16-24 UNF Thread and O-Ring



works with...

**BPA013 (-001)**


Blood Pressure Adapter with Internal 5/16-24 UNF Thread and O-Ring to 500 Series Barb, 1/8" (3.2 mm) ID Tubing

**BPA025 (-001)**


Blood Pressure Adapter with Internal 5/16-24 UNF Thread and O-Ring to 500 Series Barb, 5/32" (4.0 mm) ID Tubing

**BPA035 (-001)**

Blood Pressure Adapter with Internal 5/16-24 UNF Thread and O-Ring to 500 Series Barb, 3/16" (4.8 mm) ID Tubing

**BPA055 (-001)**

Blood Pressure Adapter with Internal 5/16-24 UNF Thread and O-Ring to 500 Series Barb, 1/4" (6.4 mm) ID Tubing



THE BPLF SERIES

CUFF CONNECTORS

Value Plastics introduces BPLF Series quarter-turn cuff connectors, expanding the range of solutions offered to blood pressure cuff manufacturers. BPLF Series products feature an ergonomic grip, facilitating gripping the fitting, connecting and disconnecting the cuff as needed. They mate perfectly with quarter-turn male o-ring connectors already found on some patient monitoring equipment.

BPLF230 (-8012)

Quarter-turn Female Blood Pressure Cuff Connector to 200 Series Barb, 1/8" (3.2 mm) Tubing ID

**BPLF240 (-8012)**

Quarter-turn Female Blood Pressure Cuff Connector to 200 Series Barb, 5/32" (4.0 mm) Tubing ID

