

### NEW 15,000 psi UHPLC microbore internal sample injectors, 1/16" Valco fittings, 0.25 mm ports (.010")

Model C74X

**SPECS**  
15,000 psi liq  
50°C max  
Stainless stator  
with inert coating  
Valcon E3 rotor

Includes stainless steel  
nuts and ferrules.

Standard electric actuator:  
110 VAC for USA  
110/230 VAC to 24 VDC power supply for  
international.  
Microelectric actuator:  
24 VDC, with 110/230 VAC to 24 VDC  
power supply.



15,000 psi

Microbore

Internal sample

1/16" 0.25 mm

#### OPTIONS

0.15 mm ports (.006")

	10 nanoliters		20 nanoliters		50 nanoliters	
	Prod No	Price	Prod No	Price	Prod No	Price
Manual	C74X-1694-.01		C74X-1694-.02		C74X-1694-.05	
With pneumatic actuator	C74X-1694-.01A		C74X-1694-.02A		C74X-1694-.05A	
With standard electric actuator	C74X-1694-.01E		C74X-1694-.02E		C74X-1694-.05E	
With microelectric actuator	C74X-1694-.01EH		C74X-1694-.02EH		C74X-1694-.05EH	
Replacement valve	C74X-1694-.01D		C74X-1694-.02D		C74X-1694-.05D	
Replacement rotor	C74-16R-.01		C74-16R-.02		C74-16R-.05	
Replacement stator	C74-1C9		C74-1C9		C74-1C9	

### NEW 10,000 psi UHPLC microbore internal sample injectors, 1/16" Valco fittings, 0.25 mm ports (.010")

Model C74H

**SPECS**  
10,000 psi liq  
50°C max  
Stainless stator  
with inert coating  
Valcon E3 rotor

Includes stainless steel  
nuts and ferrules.

Standard electric actuator:  
110 VAC for USA  
110/230 VAC to 24 VDC power supply for  
international.  
Microelectric actuator:  
24 VDC, with 110/230 VAC to 24 VDC  
power supply.



10,000 psi

Microbore

Internal sample

1/16" 0.25 mm

#### OPTIONS

0.15 mm ports (.006")

	10 nanoliters		20 nanoliters		50 nanoliters	
	Prod No	Price	Prod No	Price	Prod No	Price
Manual	C74H-1694-.01		C74H-1694-.02		C74H-1694-.05	
With pneumatic actuator	C74H-1694-.01A		C74H-1694-.02A		C74H-1694-.05A	
With standard electric actuator	C74H-1694-.01E		C74H-1694-.02E		C74H-1694-.05E	
With microelectric actuator	C74H-1694-.01EH		C74H-1694-.02EH		C74H-1694-.05EH	
Replacement valve	C74H-1694-.01D		C74H-1694-.02D		C74H-1694-.05D	
Replacement rotor	C74-16R-.01		C74-16R-.02		C74-16R-.05	
Replacement stator	C74-1C9		C74-1C9		C74-1C9	

#### MORE INFORMATION

Actuators

Air ..... page 195

Microelectric .. 188-189

Standard electric ... 193

Materials

Metals..... 254-255

Polymers ..... 256

Valve rotors..... 257

Standoff

assemblies .... 205-207



Model C74H  
(C74X similar)  
1/16" ZDV fittings

CHROMalytic TECHNOlogy Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034