

Selectors – Low Pressure


Stream selector, 1/16" Valco ZDV fittings, 0.75 mm ports (.030")

Model C25Z

- Low pressure
- Stream selector
- 10-32 ZDV
- 1/16" 0.75 mm

Includes Valco ZDV PEEK nuts and ferrules.

Standard electric actuator:
110 VAC for USA
110/230 VAC to 24 VDC power supply for international.
Microelectric actuator:
24 VDC, with 110/230 VAC to 24 VDC power supply.




SPECS

100 psi gas/ 250 psi liq
75°C max
PPS stator
Valcon E2 rotor

OPTIONS

- 4 and 12 positions available
- 2", 3", 4", and 6" standoffs
- Other polymeric materials are available. Consult the factory.

	6 Position		8 Position		10 Position		14 Position	
	Prod No	Price	Prod No	Price	Prod No	Price	Prod No	Price
Manual	C25Z-3186		C25Z-3188		C25Z-3180		C25Z-31814	
With pneumatic act.	C25Z-3186A		C25Z-3188A		C25Z-3180A		C25Z-31814A	
With std electric act.	C25Z-3186E		C25Z-3188E		C25Z-3180E		C25Z-31814E	
With microelectric act.	C25Z-3186EMH		C25Z-3188EMH		C25Z-3180EMH		C25Z-31814EMH	
Replacement valve	C25Z-3186D		C25Z-3188D		C25Z-3180D		C25Z-31814D	
Replacement rotor	C15-310		C15-310		C15-310		C25Z-325	
Replacement stator	C25Z-386		C25Z-388		C25Z-380		C25Z-38-14	



Model C25Z
10 positions
1/16" ZDV fittings

OPTIONAL FLOWPATH

Model C25ZF, the flow-through version, is similar to the C25Z but its non-selected streams continue flowing through individual outlets, instead of being dead-ended. 3, 4, 5, 6, and 7 positions are available.

Consult the factory for C25ZF prices and information.

