

Microbore centered port injector, 1/16" Valco fittings, 0.25 mm ports (.010")

Model C3

SPECS
5000 psi liq
75°C max
 Metal stator
 Valcon H rotor

5000 psi liq
50°C max
 PAEK stator
 Valcon E rotor

Includes stainless steel nuts and ferrules. Valves with PAEK stators have PEEK nuts and ferrules.

Includes one 5 µl loop of the stator material. Includes syringe fill port for 22 gauge 3/4" and 2" needle.


Standard electric actuator:
 110 VAC for USA
 110/230 VAC to 24 VDC power supply for international.

Microelectric actuator:
 24 VDC, with 110/230 VAC to 24 VDC power supply.



- 5,000 psi**
- Microbore**
- Centered port**
- 1/16"**
- 0.25 mm**

N60 stainless stator		PAEK stator	
Prod No	Price	Prod No	Price
Manual	C3-1006	C3-1346	
With pneumatic actuator	C3-1006A	C3-1346A	
With standard electric actuator	C3-1006E	C3-1346E	
With microelectric actuator	C3-1006EH	C3-1346EH	
Replacement valve	C3-1006D	C3-1346D	
Replacement rotor	C2-10R6	C2-13R6	
Replacement stator	C3-1C06	C3-1C46	



OPTIONS

- Titanium and Hastelloy stators available



Microbore vertical port injector, 1/16" Valco fittings, 0.25 mm ports (.010")

Model C2V

SPECS
5000 psi liq
75°C max
 Metal stator
 Valcon H rotor

5000 psi liq
50°C max
 PAEK stator
 Valcon E rotor

Includes stainless steel nuts and ferrules. Valves with PAEK stators have PEEK nuts and ferrules.

Includes one 5 µl loop of the stator material.


Standard electric actuator:
 110 VAC for USA
 110/230 VAC to 24 VDC power supply for international.

Microelectric actuator:
 24 VDC, with 110/230 VAC to 24 VDC power supply.



- 5,000 psi**
- Microbore**
- Vertical port**
- 1/16"**
- 0.25 mm**

N60 stainless stator		PAEK stator	
Prod No	Price	Prod No	Price
Manual	C2V-1006	C2V-1346	
With pneumatic actuator	C2V-1006A	C2V-1346A	
With standard electric actuator	C2V-1006E	C2V-1346E	
With microelectric actuator	C2V-1006EH	C2V-1346EH	
Replacement valve	C2V-1006D	C2V-1346D	
Replacement rotor	C2-10R6	C2-13R6	
Replacement stator	C2V-1C06	C2V-1C46	



OPTIONS

- Titanium and Hastelloy stators available

