

NEW Integrated motor/valve, 1/4-28 fitting details for 1/16" tubing, 0.75 mm ports (.030")

Model C62

SPECS
 100 psi gas/ 250 psi liq
 75°C max
 PPS stator
 Valcon E2 rotor

Includes multicolored Cheminert flangeless fittings for 1/16" tubing.

Sample loops are not included with valves. Order separately.

Low pressure

Integrated

1/4-28 Internal

1/16" 0.75 mm

CE



4 Port

Prod No



6 Port

Prod No



8 Port

Prod No



10 Port

Prod No

With motor, sensor, & controller
 With motor and sensor only

C62-3184I
 C62-3184I-S

C62-3186I
 C62-3186I-S

C62-3188I
 C62-3188I-S

C62-3180I
 C62-3180I-S

NEW Integrated motor/valve, 1/4-28 fitting details for 1/8" tubing, 1.50 mm ports (.060")

Model C62

SPECS
 100 psi gas/ 250 psi liq
 75°C max
 PPS stator
 Valcon E2 rotor

Includes multicolored Cheminert flangeless fittings for 1/8" tubing.

Sample loops are not included with valves. Order separately.

Low pressure

Integrated

1/4-28 Internal.

1/8" 1.50 mm

CE

4 Port

Prod No

6 Port

Prod No

8 Port

Prod No

10 Port

Prod No

With motor, sensor, & controller
 With motor and sensor only

C62-6184I
 C62-6184I-S

C62-6186I
 C62-6186I-S


C62-6188I
 C62-6188I-S

C62-6180I
 C62-6180I-S

Sample loops

for Model C62

Loops include flangeless fittings with natural color nuts. Loops less than 500 µl are made from 1/16" OD tubing; loops 500 µl or greater are made from 1/8" OD tubing.



Model C62
 1/4-28 fittings

	FEP	PTFE	PEEK
Volume	Prod No	Prod No	Prod No
20 µl	CFSL20FEP	CFSL20TF	CFSL20PK
50 µl	CFSL50FEP	CFSL50TF	CFSL50PK
100 µl	CFSL100FEP	CFSL100TF	CFSL100PK
250 µl	CFSL250FEP	CFSL250TF	CFSL250PK
500 µl	CFSL500FEP	CFSL500TF	CFSL500PK
1 ml	CFSL1KFEP	CFSL1KTF	CFSL1KPK
2 ml	CFSL2KFEP	CFSL2KTF	CFSL2KPK