

Closemount Hardware

If a valve is not going to be heated beyond the temperature range of the actuator, closemount hardware often makes the cleanest installation.

Closemount hardware

for manual valves

If you have a Valco W Type valve with no hardware and want a knob on it, or if you are converting an air or electrically actuated two position valve to manual use, this is what you need. There are two versions: one for valves with threaded mounting holes and one for valves with unthreaded mounting holes. (If your valve has no mounting holes, you will have to use it with a standoff.)

Description	Prod No
For valves with threaded mounting holes	WMMA
For valves with unthreaded mounting holes	WMMA10



Closemount hardware

for actuators

Order the appropriate closemount hardware if you want to change your valve and actuator from a standoff to a closemount connection. Two mounting screws are included. If air and standard electric actuators require different mounting screws, two of each screw are included with the closemount hardware.

Description	Prod No
-------------	---------

Air or standard electric actuators

For Valco two position valves with 1 or 2 mounting holes	CMH
For Valco two position valves with no mounting holes	CMHMP
For Valco multiposition valves	CMHMP
For Cheminert valves high pressure design	CMH11H
For Cheminert valves low pressure design (low pressure design includes required adapter)	CMH11L



Microelectric actuators

For Valco two position valves with 1 or 2 mounting holes	CMH12H
For Valco two position valves with no mounting holes	CMH12H
For Valco multiposition valves (UW and MW Types only)	CMH13
For Cheminert two position valves high pressure design	CMH12H
For Cheminert two position valves low pressure design (low pressure design includes required adapter)	CMH12L
For Cheminert multiposition valves high pressure design	CMH13H
For Cheminert multiposition valves low pressure design (low pressure design includes required adapter)	CMH13L

TECH TIP

If you need the **actuator as well as the hardware**, you can order it complete with the appropriate hardware or with the required standoff already installed.

Actuators

Air pages 194-195
 Microelectric 188-191
 Standard elec 193