



GC Capillary Columns from VICI Metronics

ValcoBond® and ValcoPLOT® capillary columns meet the highest quality standards for resolution, retention characteristics, inertness, bleed, and reproducibility.

ValcoBond® Capillary Columns

- Individually tested
- High temperature range
- Competitive pricing

We use proprietary liquid phase processing to produce low bleed characteristics while maintaining identical retention characteristics to the phases you are used to.

ValcoPLOT® Capillary Columns

- Widest polarity range
- Faster than micropacked
- Water tolerant

Now you can reduce run time by replacing your packed columns with ValcoPLOT HayeSep capillary PLOT columns, with phases available only from VICI. Our proprietary phase processing produces the first capillary PLOT columns with characteristics identical to HayeSep packed columns.

VALCOBOND PHASES

pages 235 - 238

VB-FLUORO <i>NEW!</i>	Bonded fluorosilicone phase
VB-1	100% dimethylpolysiloxane
VB-5	(5%-Phenyl)-methylpolysiloxane
VB-35	(35%-Phenyl)-methylpolysiloxane
VB-50/608	(50%-Phenyl)-methylpolysiloxane
VB-624	(6% Cyanopropyl-phenyl)-methylpolysiloxane
VB-1701	(14% Cyanopropyl-phenyl)-methylpolysiloxane
VB-Wax	Polyethylene glycol (PEG)



VALCO PLOT PHASES

pages 240 - 244

ValcoPLOT Molesieve 5Å	ValcoPLOT A	High purity Divinylbenzene/ethyleneglycoldimethacrylate
ValcoPLOT Metal Molesieve 5Å	ValcoPLOT B	Divinylbenzene/polyethyleneimine
ValcoPLOT Alumina KCl	ValcoPLOT C	Divinylbenzene/acrylonitrile
ValcoPLOT Alumina Na₂So₄	ValcoPLOT D	High purity Divinylbenzene
	ValcoPLOT N	Divinylbenzene/ethyleneglycoldimethacrylate
	ValcoPLOT P	Divinylbenzene/styrene
	ValcoPLOT Q	Divinylbenzene
	ValcoPLOT R	Divinylbenzene/N-vinyl-2-pyrrolidinone
	ValcoPLOT S	Divinylbenzene/4-vinyl-pyridine

PRODUCTS FOR GC

Other useful products for gas chromatography include:

1/32" ultra low mass external unions . . .	p. 19
FS adapter ferrules	16,17
GC detectors	222-225
GC injection valves	102-111
GC stream selectors	122-131
Gas purifiers	227
Helium and nitrogen purifiers	226
Inlet discs (injector nuts for HP 6890 and 5890)	19
Reduced breakdown injection port liners	245

TO ORDER

For prices or more information about your specific application, contact VICI Metronics:

Toll-free 877-737-1887
Tel 360-697-9199
Fax 360-697-6682

columns@vici.com

CHROMalytic TEChnology Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034