


**Microbore column end fittings**

(1.0 mm – 2.0 mm column ID)


Standard material is Type 316 stainless.

	Bore	Column ID	Without frit		Removable 2µ frit	
			Prod No	Price	Prod No	Price
<b>External column end fittings</b>						
1/16" to 1/16"	0.25 mm	1.0 mm	ECEF111.0		ECEF111.0F	
1/8" to 1/16"	0.25 mm	1.0 mm	ECEF211.0		ECEF211.0F	



**Microbore external column end fitting**  
(ECEF211.0F)

	Bore	Column ID	Without frit		Removable 2µ frit	
			Prod No	Price	Prod No	Price
<b>Internal column end fittings</b>						
1/16" to 1/32"	0.25 mm	1.0 mm	CEF1.5		CEF1.5F	
1/16" to 1/16"	0.25 mm	1.0 mm	CEF1		CEF1F	
1/8" to 1/32"	0.25 mm	1.0 mm	CEF2.51.0		CEF2.51.0F	
1/8" to 1/16"	0.25 mm	1.0 mm	CEF211.0		CEF211.0F	
1/8" to 1/16"	0.25 mm	2.0 mm	CEF212.0		CEF212.0F	



**Microbore internal column end fitting**  
(CEF1F)

**NANOBORE COLUMN END FITTINGS**

See our complete line of 100 µm and 150 µm bore fittings on page 62.

- 100 µm = .004"
- 150 µm = .006"
- 0.25 mm = .010"
- 0.50 mm = .020"
- 0.75 mm = .030"
- 1.0 mm = .040"
- 1.5 mm = .060"
- 2.0 mm = .080"
- 4.6 mm = .180"
- 6.0 mm = .236"
- 6.4 mm = .253"
- 7.0 mm = .275"
- 10.0 mm = .400"
- 27.0 mm = 1.08"
- 1/32" = 0.8 mm
- 1/16" = 1.6 mm
- 1/8" = 3.2 mm
- 1/4" = 6.4 mm
- 3/8" = 9.5 mm
- 1/2" = 12.7 mm