



Open end wrenches

Size	For use with	Prod No	Price
3/16" x 1/4"	1/32" and 1/16" nuts	OEW	
3/8" x 7/16"	1/8" nuts	OEW-2	
1/2" x 9/16"	1/4" nuts	OEW-3	



Pin vise and drill index


The drill index has drills sized from 0.0135" to 0.039" (0.34 to 1 mm). These are useful tools when a fused silica tube breaks in a union (see Tech Tip on the facing page) and for enlarging the inner diameter of fused silica adapters.

Prod No	Price
PV	

Template

This tool is just what you need when you're working out plumbing and valve switching schematics. It features templates for two position valves with 4, 6, 8, and 10 ports with indications of both positions, as well as various flow symbols. For added convenience, the sides are edged with metric and inch rulers.

Prod No	Price
TEMPLATE1	



MORE INFORMATION

- Tools for valves
- Pencil magnet p 210
- Valve spanner handle 211
- Tightening tools for PEEK fittings 67
- Tubing accessories 90