

### Mixing tees

1/4-28



Includes flangeless 1/4-28 fittings for tubing OD indicated.

Tubing		PEEK		CTFE	
OD	Bore	Prod No	Price	Prod No	Price
1/16"	0.75 mm	CM1XPK		CM1XKF	
1/8"	1.50 mm	CM2XPK		CM2XKF	




### Glass connectors

1/4-28 female to 1/4" glass tube

Glass connectors join a Cheminert tube end fitting to 1/4" OD glass tubing. They are available as individual components or as complete assemblies. Assemblies include a bushing and nut, a polypropylene or CTFE split ring, and a 1/4" OD x 3-1/4" long piece of 1 mm or 2 mm ID glass tube. This connector works only with our glass tubes.

Description	Acetal		CTFE	
	Prod No	Price	Prod No	Price
<b>Complete assembly</b>				
1 mm ID glass tubing	CGC41		CGC41KF	
2 mm ID glass tubing	CGC42		CGC42KF	
<b>Components</b>				
Bushing	CGCB		CGCBKF	
Nut	CGCN		CGCNKF	
Glass rod, 3-1/4" long				
1 mm ID	CGCG41		-	
2 mm ID	CGCG42		-	
Split rings (package of 5)	CGCR		CGCRKF	



Glass connector