

Pipe adapters

1/4-28 to NPT

Versions adapt male or female 1/4-28 fittings to male or female NPT.

| NPT | Bore | PEEK | | CTFE | |
|------------------------------------|--------|---------|-------|---------|-------|
| | | Prod No | Price | Prod No | Price |
| Female 1/4-28 to male NPT | | | | | |
| 1/8" | 1.5 mm | CPA2PK | | CPA2KF | \$18 |
| 1/4" | 1.5 mm | CPA4PK | | CPA4KF | |
| Male 1/4-28 to male NPT | | | | | |
| 1/8" | 1.5 mm | CEPA2PK | | CEPA2KF | |
| 1/4" | 1.5 mm | CEPA4PK | | CEPA4KF | |
| Female 1/4-28 to female NPT | | | | | |
| 1/8" | 1.5 mm | CFPA2PK | | CFPA2KF | |
| 1/4" | 1.5 mm | CFPA4PK | | CFPA4KF | |



Cheminert 1/4-28 to Valco 10-32 ZDV adapter


This adapter permits Valco 10-32 fittings to be installed into any 1/4-28 fitting detail. (Nut and ferrule are not included.)

| Description | Bore | Prod No | Price |
|--------------|---------|---------|-------|
| Port adapter | 0.50 mm | ZLCA1PK | |

One-piece fingertight column coupler

Choose from a variety of coupler IDs, indicated by the color of the sleeve (which parallels the color-coding of our PEEK tubing on page 89). A unique feature of this column coupler is that it adapts automatically to fit all pilot lengths – Valco, Waters, Upchurch, Rheodyne, etc. Since the tubing bottoms out in any fitting detail, added void volume is minimal. Material is PEEK.

| Color | Bore | Prod No | Price |
|--------|------------|----------|-------|
| Red | 0.13 mm ID | JR-26501 | \$20 |
| Yellow | 0.17 mm ID | JR-26502 | 20 |
| Blue | 0.25 mm ID | JR-26503 | 20 |
| Orange | 0.50 mm ID | JR-26504 | 20 |



| |
|-----------------|
| 0.25 mm = .010" |
| 0.50 mm = .020" |
| 0.75 mm = .030" |
| 1.0 mm = .040" |
| 1.5 mm = .060" |
| 2.0 mm = .080" |
| 4.6 mm = .180" |
| 6.0 mm = .236" |
| 6.4 mm = .253" |
| 7.0 mm = .275" |
| 10.0 mm = .400" |
| 27.0 mm = 1.08" |
| 1/32" = 0.8 mm |
| 1/16" = 1.6 mm |
| 1/8" = 3.2 mm |
| 1/4" = 6.4 mm |
| 3/8" = 9.5 mm |
| 1/2" = 12.7 mm |