

LARGE BORE



MALE
PG-27




FEMALE
PG-27




PLUGS AND CAPS
PG-27




MLB035	NUMBER	MATERIAL
Male Large Bore Connector to 500 Series Barb, 3/16" (4.8 mm) ID Tubing	MLB035-1	White Nylon
	MLB035-8012	Medical Grade White ABS
	MLB035-9	Clear Polycarbonate




FLB035	NUMBER	MATERIAL
Female Large Bore Connector to 500 Series Barb, 3/16" (4.8 mm) ID Tubing	FLB035-1	White Nylon
	FLB035-8012	Medical Grade White ABS
	FLB035-9	Clear Polycarbonate




MLB055	NUMBER	MATERIAL
Male Large Bore Connector to 500 Series Barb, 1/4" (6.4 mm) ID Tubing	MLB055-1	White Nylon
	MLB055-8012	Medical Grade White ABS
	MLB055-9	Clear Polycarbonate




FLB055	NUMBER	MATERIAL
Female Large Bore Connector to 500 Series Barb, 1/4" (6.4 mm) ID Tubing	FLB055-1	White Nylon
	FLB055-8012	Medical Grade White ABS
	FLB055-9	Clear Polycarbonate




MLB065	NUMBER	MATERIAL
Male Large Bore Connector to 500 Series Barb, 5/16" (8.0 mm) ID Tubing	MLB065-1	White Nylon
	MLB065-8012	Medical Grade White ABS
	MLB065-9	Clear Polycarbonate



FLB065	NUMBER	MATERIAL
Female Large Bore Connector to 500 Series Barb, 5/16" (8.0 mm) ID Tubing	FLB065-1	White Nylon
	FLB065-8012	Medical Grade White ABS
	FLB065-9	Clear Polycarbonate




MLBC7	NUMBER	MATERIAL
Male Large Bore Connector with Integral Lock Ring to Male Slip Luer (May be Used with Separate Rotating Lock Ring; FSLR)	MLBC7-9002	Rad. Stable Polycarbonate



LARGE BORE CONNECTORS:
Plugs and Caps

MLBC8	NUMBER	MATERIAL
Male Large Bore Connector with Integral Lock Ring to Female Thread-Style Locking Luer	MLBC8-9002	Rad. Stable Polycarbonate



FLBP	NUMBER	MATERIAL
Female Large Bore Cap	FLBP-1	White Nylon
	FLBP-8012	Medical Grade White ABS
	FLBP-9	Clear Polycarbonate



MLBP	NUMBER	MATERIAL
Male Large Bore Plug	MLBP-1	White Nylon
	MLBP-8012	Medical Grade White ABS
	MLBP-9	Clear Polycarbonate



LARGE BORE

HRomalytic www.chromtech.net.au sales@chromtech.net.au
 ABN 14 643 445 058 **PTY LTD** Tel: (03) 9762 2024 Fax: +61 3 9761 1169
TECH *nology*
 Australian Distributors: **TECH** *nology* ... for ALL Your chromatography supplies!