

BLOOD PRESSURE



BPF SERIES PG-27



BAYONET MALE CUFF CONNECTORS & DINAMAP® MONITOR ADAPTERS PG-28



BPLF 1/4 TURN CUFF CONNECTORS PG-28



CUFF CONNECTORS PG-28



THE **BPF** SERIES

From Value Plastics, the reference standard for healthcare fluid management products: The BPF Series bayonet-style female quick connector. BPF series quick connect fittings are designed to withstand abuse, provide economy and ease of use over other alternatives, and comply with ISO Standard BS EN 1060-3.

Ergonomic Shape

- Easy to hold, even with gloves
- Won't irritate patients in case of incidental contact

Large Guarded Button

- Easy to use, even with gloves
- Easy to understand, convenient operation
- Protected against accidental disconnect

Push To Connect

- Easiest connector to use, even with gloves

O-Ring Seal

- Never gets more difficult to put together
- Consistent performance throughout product life

Bayonet Rotates In Connector

- Prevents kinked tubes
- Facilitates easy connections to patient cuff

Unique Spring-Loaded Bayonet Retainer

- Clicks when assembled to signal secure connection
- Resists accidental disconnects

Rugged Polymer Construction

- Withstands tens of thousands of cycles
- Withstands hospital sanitizing procedures

Unimpeded Flow Path

- Better acoustics for accurate readings