

# TAPERED SEAL




**SMALL BORE  
PG-30**



**LARGE BORE MALE  
PG-31**



**LARGE BORE FEMALE  
PG-31**



**PLUGS AND CAPS  
PG-31**



## **SMALL BORE CONNECTORS: Plugs and Caps**

### **TSC2M013 (-6)**

Male Taper Seal Non-Luer Connector to 500 Series Barb, 1/8" (3.2 mm) ID Tubing



### **TSC2M035 (-6)**

Male Taper Seal Non-Luer Connector to 500 Series Barb, 3/16" (4.8 mm) ID Tubing



### **TSC2F007 (-6)**

Female Taper Seal Non-Luer Connector to 500 Series Barb, 3/32" (2.4 mm) ID Tubing



### **TSC2F013 (-6)**

Female Taper Seal Non-Luer Connector to 500 Series Barb, 1/8" (3.2 mm) ID Tubing



### **TSC2F035 (-6)**

Female Taper Seal Non-Luer Connector to 500 Series Barb, 3/16" (4.8 mm) ID Tubing



### **TSC2FB230 (-6)**

Female Taper Seal Non-Luer with 1/4-28 UNF Panel Mount to 200 Series Barb, 1/8" (3.2 mm) ID Tubing



### **TSC2MP (-6)**

Male Taper Seal (Non-Luer) Plug with Tether Recess



**LARGE BORE CONNECTORS:**

*Male*

**LARGE BORE CONNECTORS:**

*Female*

**Material Suffix Key**

- 1 White Nylon
- 4 Green Nylon
- 6 Natural Polypropylene
- 9 Clear Polycarbonate


- 10 Acetal
- 5 Blue Nylon
- 6005 Animal-Free Polypropylene
- 9002 Radiation Stable Polycarbonate

- 2 Black Nylon
- 06 Orange Nylon
- 8003 Natural ABS
- 11A Natural Kynar PVDF

- 3 Red Nylon
- 07 Yellow Nylon
- 81 Lace White ABS


**MLB035 (-1 -8012 -9)**

Male Large Bore Connector to 500 Series Barb, 3/16" (4.8 mm) ID Tubing




**FLB035 (-1 -8012 -9)**

Female Large Bore Connector to 500 Series Barb, 3/16" (4.8 mm) ID Tubing




**MLB055 (-1 -8012 -9)**

Male Large Bore Connector to 500 Series Barb, 1/4" (6.4 mm) ID Tubing




**FLB055 (-1 -8012 -9)**

Female Large Bore Connector to 500 Series Barb, 1/4" (6.4 mm) ID Tubing




**MLB065 (-1 -8012 -9)**

Male Large Bore Connector to 500 Series Barb, 5/16" (8.0 mm) ID Tubing




**FLB065 (-1 -8012 -9)**

Female Large Bore Connector to 500 Series Barb, 5/16" (8.0 mm) ID Tubing



**MLBC7 (-9002)**

Male Large Bore Connector with Integral Lock Ring to Male Slip Luer (May be Used with Separate Rotating Lock Ring: FSLLR)



**LARGE BORE CONNECTORS:**

*Plugs and Caps*


**FLBP (-1 -8012 -9)**

Female Large Bore Cap



**MLBC8 (-9002)**

Male Large Bore Connector with Integral Lock Ring to Female Thread-Style Locking Luer



**MLBP (-1 -8012 -9 -9002)**

Male Large Bore Plug



TAPERED SEAL



our online ordering for component parts, now available 24 hours a day, 7 days a week, with no \$ minimums! You can find every catalog part at:

[www.valueplastics.com](http://www.valueplastics.com)

(Online ordering is not available in countries served by our regional partners)