


HPLC Detector Lamps and UV-VIS



DON't get "trapped" by high OEM labor prices to replace HPLC **Detector Lamps** and **Pump parts**
 It IS relatively simple to do ! **Do-It-Yourself !**
 NO Harder than replacing a car head-light or a globe in a light socket
 WARRANTIES **Don't COVER** these consumables **ANYWAY!**

SPECIAL
\$1050 each !
Agilent/HP 1100
Diode Array Detector

Chromalytic OFFER other OEM Replacement Lamps AND at reasonable prices TOO! OFTEN 2/3 of OEM part prices

HPLC Detector Lamps

We can supply; **Deuterium (D2), Mercury (Hg), Tungsten (W), Xenon (Xe), Zinc (Zn) and Cadmium (Cd)** sources.
 Some lamps are pre-aligned.
 However some lamps (with holders) have to be returned to our facility for retrofitting.
 PLEASE Enquire - before ordering

HPLC Pump Spares Piston seals

- most are made from Graphite filled PTFE jackets (black colour) and are fitted with a stainless steel spring that acts as the seal energiser.

However, recent introductions have increased seal life considerably in specific solvent conditions. These jackets are made either of UPC (Ultra high molecular weight polyethylene) which appears as a clear (white) colour, or a special proprietary material that is gold (yellow)coloured.

Pistons (or plungers) are usually manufactured from Sapphire, but some Ruby pistons are available, which increases usage by 10%.

Check Valves and Check Valve Cartridges - Enquire !

Specifications, materials, and suppliers may change at any time ask for a quotation and enquire on availability - some models are NOW obsolete

- [ACS, Biorad, Dupont, Eldex, Jasco](#)
- [Cecil Dionex, Knauer, LKB-Parmacia, Merck-Hitachi](#)
- [Beckman](#)
- [Waters](#)
- [LDC, Milton Roy \(TSP\)](#)
- [Gilson/Rainin](#)
- [Kratos/ ABI](#)
- [Severn Analytical](#)
- [Spectra Physics \(TSP\)](#)

100's of different model numbers for different OEM instruments

- [Agilent / Hewlett Packard](#)
- [Kontron](#)
- [SSI/ Bischoff](#)
- [Perkin Elmer](#)
- [BAS, Erma, Micromeritics, Shodex](#)
- [ATI Unicam \(Pye Unicam Philips\)](#)
- [Shimadzu](#)
- [Varian](#)

Also lamps for UV-VIS and AA

This is ONLY a part list of those available - Enquire before ordering !
 Enquire about lamps and pump parts for other OEMs and for the various model numbers