

# SRI Model 310MM Cannabis and Hemp Gas Chromatograph with Built-in Hydrogen Generator and Air Compressor

**The SRI Model 310MM Gas Chromatograph (GC)** ( part# 0310-0091 ) is the ideal GC for many labs and even individuals seeking to test cannabis and hemp.

**Measure** cannabinoids, terpenes, and residual solvents in cannabis concentrates.

**Low Cost.** At \$11,134.00, the Model 310MM is a fraction of the cost of other gas chromatographs and analytical instruments. *(2022 pricing, prices subject to change, consult most recent price list.)*

**Portable.** You can fly with the GC as airline baggage.

**Rugged.** The Model 310MM is designed for rugged use in the lab or the field.

**Fast.** Ten minutes or less for a cannabinoid analysis. Sixteen minutes for terpenes and cannabinoids in a single test.

**No gas cylinders required.** The built-in Hydrogen generator and air compressor provide all required gases, running only on distilled water.

**Powerful.** Two columns let you separate gases or liquids without changing any hardware while the sensitive FID detector detects down to .01% THC.

**Easy.** Anybody can learn to operate the GC in a few hours. **You do not have to be a scientist.**

**Free Software Included.** And you can install the software on as many computers as you like. Updates are free.

**Two Year Warranty.** One year on H2 cell.

**Free Training** and United States based technical support.

**Videos** Watch our many YouTube videos. Click the icon at the bottom of the srigc.com webpage.

**To Order:**

SRI Model 310MM GC

Part# 0310-0091

