



Zhejiang ALWSCI Technologies Co., Ltd.

Yuexing Technology Innovation Park, No. 11 Fanrong Road, Binhai New City,
Shaoxing 312366, Zhejiang P.R., China.

Product Certificate

CATALOG NUMBER C0000342 LOT NUMBER 00E0923VJ

GC/MS Agilent 7890/5975

TOC Analyzer EST Vocarb 3000

PROCESS DATE 01/23/21

PARAMETERS

DESCRIPTION 40mL 27.5x95mm Clear Glass EPA/TOC Vial 24-400 White Open Top PP Screw
Cap with 22mm Natural PTFE/White Silicone 3.0mm thick Septa (EPA Quality)
and Dust Cover.

This is to certify that the item(s) as stated above have been processed on Agilent 7890/5975 and EST Vocarb 3000

ANALYSIS

Total vial contribution after blank subtraction is as listed:

ANALYTE	SPECIFICATION	RESULT
Total Organic Carbon	<10ug/L	<10ug/L

ALWSCI Technologies certifies that the item(s) in above referenced shipment were processed, tested, and meets the specifications above for this product.

Approved By: _____

Mary Ma 01/23/2021

Quality Assurance