

Hub-Cap Mobile Phase Accessories

Simplify Mobile Phase Delivery with the Hub-Cap Filter Kit

Introducing our new Hub Cap filter kit! The Hub-Cap filter allows you to simultaneously transfer and filter your mobile phases. The bottle tops and adaptors are designed to fit securely on 4-liter solvent bottles and eliminate messy, loose-fitting parafilm or foil wraps. Tidy up your mobile phase delivery—try a Hub-Cap today!



Keep your mobile phase lines under control—use Hub-Cap bottle tops instead of parafilm, aluminum foil, or tape on your mobile phase reservoirs.

Hub-Cap Filter Kit

new!

Transfer and filter mobile phases in one easy step.

Kit includes: bottle adapter, bottle adapter nut, filter inlet cap, grid support, vacuum hose barb, tube compression fitting, 47mm grid, 47mm .22 μ m filter membrane, 47mm .45 μ m filter membrane, 1/4" OD x 1/8" ID ultra chemical resistant, Teflon® FEP lined Tygon® tubing (3'), 6" x 6" box with shrink wrap insert



cat. #26395

| Description | qty. | cat.# |
|--|---------|-------|
| Hub-Cap Filter Kit | kit | 26395 |
| Replacement Membrane Filters | qty. | cat.# |
| Polypropylene Membrane Filters, 47mm, 0.45 μ m | 100-pk. | 26396 |
| Polypropylene Membrane Filters, 47mm, 0.22 μ m | 100-pk. | 26397 |
| Nylon Membrane Filters, 47mm, 0.45 μ m | 100-pk. | 26398 |
| Nylon Membrane Filters, 47mm, 0.22 μ m | 100-pk. | 26399 |

Hub-Cap 4 Liter Bottle Tops

Most bottles use a GL45 cap. New Hub-Cap bottle tops are a great way to neatly keep your mobile phase lines where they belong. Use instead of parafilm, aluminum foil, or tape on your mobile phase reservoirs.



cat. #26541

| Description | qty. | cat.# |
|---|-------|-------|
| Hub-Cap (assembly of the bottle cap and plug) | kit | 26541 |
| Hub-Cap Multi-pack | 3-pk. | 26542 |

Hub-Cap Adapters

Allow the use of the Opti-Cap™ with 4-liter solvent bottles.



cat. #26538

cat. #2654

| Description | qty. | cat.# |
|-------------------------------|-------|-------|
| Hub-Cap Adapter | ea. | 26538 |
| Hub-Cap Adapter Multi-pack | 3-pk. | 26539 |
| Hub-Cap Adapter and Opti-Cap™ | kit | 26540 |