

Now Available! Capillary Stainless Steel Tubing Assemblies

For Agilent HPLC systems

By Becky Wittrig, Ph.D, HPLC Product Marketing Manager



26525
Capillary SS Tubing
130mm x 0.17mm ID
with fittings

- Precut, micropolished tubing and preassembled fittings for quick, easy maintenance of your systems.
- Designed and tested for Agilent HPLC systems.
- Restek offers a full range of high-quality replacement parts for your HPLC systems.

Capillary Stainless Steel Tubing Assemblies for Agilent HPLC Systems

| Description | Similar to Agilent part # | qty. | cat.# |
|---|---------------------------|--------|-------|
| Capillary SS Tubing, 130mm x 0.17mm ID, with fittings | 01090-87305 | ea. | 26525 |
| Capillary SS Tubing, 800mm x 0.17mm ID, with fittings | 01078-87305 | ea. | 26526 |
| Capillary SS Tubing, 180mm x 0.17mm ID, with fittings | G1313-87305 | ea. | 26527 |
| Capillary SS Tubing, 700mm x 0.25mm ID, with fittings | 01018-67305 | ea. | 26528 |
| Capillary SS Tubing, 700mm x 0.25mm ID, with fittings | 01078-87306 | ea. | 26529 |
| Seat Capillary, SS Tubing, 0.17mm ID | 01078-87303 | ea. | 26530 |
| Capillary SS Tubing, 105mm x 0.17mm ID | 5021-1816 | ea. | 26531 |
| Mixing Capillary Assembly | G1312-67302 | ea. | 26532 |
| Capillary SS Tubing, Valve to Metering Head | G1313-87301 | ea. | 26533 |
| Capillary SS Tubing, 150mm x 0.17mm ID | 5021-1817 | ea. | 26534 |
| Capillary SS Tubing, 280mm x 0.17mm ID | 5021-1818 | ea. | 26535 |
| Capillary SS Tubing, 400mm x 0.17mm ID | 5021-1819 | ea. | 26536 |
| 1/16" SS Fitting, Front and Back Ferrules | 5062-2418 | 10-pk. | 26537 |



26530
Seat Capillary, SS Tubing
0.17mm ID



26527
Capillary SS Tubing
180mm x 0.17mm ID
with fittings



26533
Capillary SS Tubing
Valve to Metering Head



26529
Capillary SS Tubing
700mm x 0.25mm ID
with fittings



26532
Mixing Capillary
Assembly



26537
1/16" SS Fitting
Front and Back Ferrules

Get More!

For more information on
HPLC Accessories,
visit us online at
www.restek.com/HPLC



HROMalytic Chromatography Products
Australian Distributors **ECH**nology '07
www.chromtech.net.au E-mail: info@chromtech.net.au Tel: +61 3 9762 2034 Fax: +61 3 9761 1169