

Hub-Cap Bottle Tops and Adapters

Mobile Phase Accessories

- Specifically designed for 4-liter bottles!
- Easiest, cleanest way to helium sparge and deliver mobile phases.
- Use Hub-Cap bottle tops instead of parafilm, aluminum foil or tape.



HRMalytic **Chromatography**
Products '08
Australian Distributors **ECH**nology
www.chromtech.net.au E-mail : info@chromtech.net.au Tel : +61 3 9762 2034 Fax : +61 3 9761 1169

RESTEK

RESTEK

Chromatography Products

www.restek.com 800-356-1688 • 814-353-1300

Hub-Cap Bottle Tops and Adapters

Hub-Cap 4 Liter Bottle Tops

Description	qty.	cat.#
Hub-Cap (assembly of the bottle cap and plug)	kit	26541
Hub-Cap Multi-pack	3-pk.	26542



Hub-Cap Adapters


Allows the use of the Opti-Cap™ with 4-liter solvent bottles.

Description	qty.	cat.#
Hub-Cap Adapter	ea.	26538
Hub-Cap Adapter Multi-pack	3-pk.	26539
Hub-Cap Adapter and Opti-Cap™	kit	26540



cat. # 26538

cat. # 26540



Use Hub-Cap bottle tops instead of parafilm, aluminum foil or tape on your mobile phase reservoirs.

Opti-Cap™ Bottle Top

The most economical way to helium-sparge and deliver HPLC mobile phases. The Opti-Cap™ top fits all standard GL-45 bottles and has two 1/8-inch holes and one 1/16-inch hole for tubing.

Description	qty.	cat.#
Opti-Cap™ (Cap and PEEK® Plug)	ea.	25300
Opti-Cap™ Kit (Opti-Cap™, 3 meters of tubing, sparging filters)	kit	25301
Opti-Cap™ Kit with 1L Bottle	kit	25302
Opti-Cap™ Kit with 2L Bottle	kit	25303
Related items and replacement parts	qty.	cat.#
Mobile Phase Sparge Filter: 2µm, stainless steel	ea.	25311
Mobile Phase Inlet Filter: 10µm	ea.	25312
Teflon® Tubing, 1/8" OD x 0.094" ID x 3m (2.4mm ID)	3m	25307
Teflon® Tubing, 1/8" OD x 0.063" ID x 3m (1.6mm ID)	3m	25306
PEEK® Plug, 1/4"-28 threads	3-pk.	25319
1L Graduated Safety-Coated Bottle – GL-45 threads	ea.	25304
2L Graduated Safety-Coated Bottle – GL-45 threads	ea.	25305



Opti-Cap™ Kit with bottle

Restek Trademarks:
Restek logo.

Other Trademarks:
Opti-Cap (Jour Research), PEEK (Vicatex plc), Teflon (E.I. du Pont de Nemours & Co., Inc.)



Lit. Cat.# 580144

© 2007 Restek Corporation.

