



## Advantage™ Series of 96-Well Collection Plates



### Advantage™ Standard Series: 1mL 96-Well Collection Plates with Round Well & Rounded Bottom

- High quality, solvent resistant polypropylene
- Industry standard, one piece 12 X 8 well design
- Individual wells have distinct alpha-numeric designation
- Compatible with all robotic systems



### Advantage™ Clear-Vu™ Series: 1mL 96-Well Collection Plates with Round Well & Rounded Bottom

- See-thru polypropylene
- High quality, solvent resistant polypropylene
- Industry standard, one piece 12 X 8 well design
- Individual wells have distinct alpha-numeric designation
- Compatible with all robotic systems

### Advantage™ Small Volume Collection Plate with Round well

