



Stock Rubber Stoppers



We hold a large range of rubber stoppers available for immediate dispatch from our warehouse

INJECTION STOPPERS / INFUSION STOPPERS

Adelphi Part No.	Description	West Art.	Formulation	Final treatment	Lamination	Silicone	Packaging	Pack Qty	Carton Qty	MOQ
INJ13TB3WRSNP1	13mm Injection Stopper FluroTec®	1358	4023/50 Grey	NovaPure® RS	Daikyo fluoro resin D	B2-40	5 x STERILizableBAG™	1,000	5,000	1,000
INJ20TB3WRSNP	20mm Injection Stopper FluroTec®	1343	4023/50 Grey	NovaPure® RS	Daikyo fluoro resin D	B2-40	2 x STERILizableBAG™	3,000	6,000	3,000
INJ13TB3WRSRS	13mm Injection Stopper FluroTec®	1358	4023/50 Grey	Westar Select® RS	Daikyo fluoro resin D	B2-40	STERILizableBAG™	2,000	10,000	2,000
INJ20TB3WRSRS	20mm Injection Stopper FluroTec®	1343	4023/50 Grey	Westar Select® RS	Daikyo fluoro resin D	B2-40	STERILizableBAG™	2,000	6,000	2,000
INJ13TBWRS	13mm Injection Stopper Flurotec®	1258	4110/40 Grey	Westar® RS	Flurotec® +	B2-40	STERILizableBAG™	5,000	20,000	Please Contact Us
INJ20TBWRS	20mm Injection Stopper Flurotec®	1243	4110/40 Grey	Westar® RS	Flurotec® +	B2-40	STERILizableBAG™	1,000	5,000	1,000
INJ13TB3NPRU	13mm FluroTec® Injection Stopper Bromobutyl	1358	4432/50 Grey	NovaPure® RU	Daikyo fluoro resin D	B2-40	Tripple Bagged RU	1,000	5,000	1,000
INJ20TB3NPRU	20mm FluroTec® Injection Stopper Bromobutyl	1343	4023/50 Grey	NovaPure® RU	Daikyo fluoro resin D	B2-40	Tripple Bagged RU	500	2,500	500
1970-0004	13mm Injection Stopper FluroTec®	S2-F451	4432/50 Grey	Westar® RS	Daikyo fluoro resin D	B2-40	STERILizableBAG™	1,000	5,000	1,000
INJW13RtS	13mm Injection Stopper	1104	4023/50 Grey	RTS™	none	A	STERILizableBAG™	4,000	16,000	4,000
INJW20RtS	20mm Injection Stopper	1071	4023/50 Grey	RTS™	none	A	STERILizableBAG™	1,500	4,500	1,500
INJW20S	20mm Injection Stopper Silicone	1185	8763/40	Standard	none	none	Foilbag	1,500	4,500	1,500
INJW32G3RTS	32mm Infusion Stopper	4440	4023/50 Grey	RTS™	none	A	STERILizableBAG™	500	1,500	500

FREEZE DRY STOPPERS

FD13TT3WRS	13mm Freeze Dry Stopper FluroTec®	1356	4023/50 Grey	Westar® RS	Daikyo fluoro resin D	none	STERILizableBAG™	2,000	10,000	Please Contact Us
FD13TTB3WRS	13mm Freeze Dry Stopper FluroTec®	1356	4023/50 Grey	Westar® RS	Daikyo fluoro resin D	B2-TR	STERILizableBAG™	2,000	10,000	2,000
FD20TT3WRS	20mm Freeze Dry Stopper FluroTec®	1346	4023/50 Grey	Westar® RS	Daikyo fluoro resin D	none	STERILizableBAG™	1,000	4,000	Please Contact Us
FD20TTB3WRS	20mm Freeze Dry Stopper FluroTec®	1346	4023/50 Grey	Westar® RS	Daikyo fluoro resin D	B2-TR	STERILizableBAG™	2,000	4,000	2,000
FD13TBrtS	13mm Freeze Dry Stopper Flurotec®	V2-F8	D713 Grey	RTS™	Daikyo fluoro resin D	B2-40	STERILizableBAG™	5,000	10,000	5,000
FD20TBrtS	20mm Freeze Dry Stopper Flurotec®	V10-F22	D713 Grey	RTS™	Daikyo fluoro resin D	B2-40	STERILizableBAG™	1,750	3,500	1,750
FDW13WRS	13mm Freeze Dry Stopper	1079	4023/50 Grey	Westar® RS	none	4.5	STERILizableBAG™	3,000	12,000	3,000
FDW20RtS	20mm Freeze Dry Stopper	1319	4023/50 Grey	EPP (RS)	none	A	STERILizableBAG™	2,000	4,000	2,000

DIAGNOSTIC STOPPERS

FDIA14WGBRTS	14mm Freeze Dry Stopper	1067	W1816 Grey	RTS™	none	B	STERILizableBAG™	3,000	9,000	3,000
FDIA18W/B	18mm Freeze Dry Stopper	1008	PH701/40 Red	Standard	none	B	Polybag	7,000	7,000	500
FDIA18WG/B	18mm Freeze Dry Stopper	1008	W1816 Grey	Standard	none	B	Polybag	7,000	7,000	500
FDIA22	22mm Freeze Dry Stopper	1010	PH701/40 Red	Standard	none	A	Polybag	3,000	3,000	500
SDIA18RtS	18mm Plug Stopper	1007	PH701/40 Red	RTS™	none	A	STERILizableBAG™	4,000	8,000	4,000

FluroTec® (FluroTec® + and Daikyo fluoro resin D)

An ETFE lamination applied to the stoppers to significantly reduce the risk of interaction of the drug product with the closure

NovaPure®

West's tightest particulate specification for stoppers, NovaPure® components have been designed and manufactured under Quality by Design principles

Ready to Sterilize (RTS™) Enhanced Performance Product Ready-to-Sterilize (EPP (RS)) and Westar® RS (WRS)

Final rinse in WFI and packed in class 100 (ISO 5) clean room into a STERILizableBAG

B2-40

Cross linked silicone coating applied to the top surface of the rubber stopper

Adelphi Healthcare Packaging
Olympus House, Mill Green Road
Haywards Heath, West Sussex
RH16 1XQ, UK

Tel: +44 (0)1444 472300
Email: sales@adelphi-hp.com
Web: www.adelphi-hp.com