

## *Repeatable performance for repeatable results*

### Features

- Precision adjustment collar with repeatable settings
- Uses interchangeable heads for crimping and decapping
- Use for cap styles such as single piece, "flip" types, spray pumps, as well as the tougher magnetic alloy caps, and many others
- Head sizes range from 7.5mm to 43mm
- Solid construction with CNC precision throughout
- Levered cam design allows for easy operation
- Comfort grip levers
- Adjustable height and lever position for optimum ergonomics

### Benefits

- Over 50 head styles kept in stock
- One unit can be used with over 100 different sizes and styles of caps
- Quick and easy set-up and change over when working with multiple sizes and styles of caps
- Recycling facilitated by separating cap from glass
- Durable enough for heavy use
- Dependable performance
- Repeatable settings
- Simple operation
- Low cost





# MICRO-XL 4300 Adjustable Bench Crimper

**PRECISE** - Crimp setting collar allows the operator to make precise incremental adjustments as fine as 0.15mm (.006") to compensate for seal or septa thickness in order to achieve an exact seal or crimp.

**REPEATABILITY** - Settings can be kept track of and quickly returned to when working with different applications. Keeping track of settings allows the operator to return to the exact setting adjustment for each size and type of seal.

**VERSATILE** - A wide variety of interchangeable heads are available to choose from to crimp or decap over 100 different sizes and styles of caps. It can also be used to crimp spray pump caps and has no problem crimping the tougher magnetic alloy caps. Heads are quickly and easily changed without the use of tools since they simply screw in and out.



**EASY TO USE** - The levered cam design operation requires very little user effort, and the "comfort-grip" lever handle can be adjusted 360 degrees to accommodate the operator.



**ERGONOMIC** - Unit height and position are adjustable to accommodate a variety of vial and bottle sizes allowing the operator to work in a comfortable position.

**QUALITY** - Unit is manufactured using stainless steel, and plated parts providing excellent wearability, and superior corrosion resistance. These parts are CNC machined from solid stock to exact tolerances to achieve repeated results with dependable performance.

**STURDY** - The polished marble base provides stability and a flat surface to work from. An optional stainless steel base is also available.



To place an order or to obtain more information, call, fax, or visit our website.

We ship all products UPS.

Billing is available with an open account, C.O.D. UPS, or Visa/Mastercard.



**KEBBY INDUSTRIES INC.**

