



## Canister Air Sampling Timer



- Program up to 12 timed events!
- Capable of both manual and automated operation.
- Perfect for either grab or time-integrated sampling.
- Long battery life; recharges conveniently using the USB port on any PC.
- All stainless steel sample flow path ensures inertness, improving accuracy.

These timers are designed to simplify both automated and manual air sampling. The easy-to-use keypad and graphic display facilitate the programming of up to 12 timed events. They offer the convenience of remote start/stop sampling and permit intermittent sampling throughout a test period. The LCD remains in sleep mode when not in use, greatly extending battery life. Timers are compatible with any canister and flow controller.

Features include: solenoid valve for sampling control, 1/4" Swagelok® inlet and outlet fittings, highly inert stainless steel flow path, and water-proof exterior for outdoor use.

Description	qty.	Cat.#
Canister Air Sampling Timer	ea.	24267

We accept Visa, MasterCard and American Express

