





**Chromalytic Amber Vials**  
 (from left to right)

CatNo	Description	Liquid in vial	Slab Size	Slab Price	Each Price
CTSCW-100V-20T-52D94H	100ml Amber glass type (2) Vial Sulphur Teated and Washed 20mm Crimp Finish	100ml	66	\$58.08	\$1.05
CTSA-23V-20T-29D60H	23ml Amber Glass Schott Vial 20mm Crimp Finish	20ml	113	94.92	\$0.88
CTSA-13V-20T-24D49H	13ml Amber Glass Schott Vial 20mm Crimp Finish	10ml	155	\$102.30	\$0.77

**Other Amber Vials Not Necessarily Stocked**  
 (refer to clear vials for similar vial image)

CatNo	Description	Liquid in vial	Slab Size	Slab Price	Each Price
CTSAW-50V-20T-43D73H	50ml Amber glass type (2)Vial Sulphur Teated and Washed 20mm Crimp Finish	50ml	112	Call	Call
CTSA-7V-20T-22D38H	7ml Amber Glass Schott Vial 20mm Crimp Finish	7ml	241	Call	Call
CTSA-2V-13T-15D32H	2ml Amber Glass Schott Vial 13mm Crimp Finish	2ml	371	Call	Call



## Hand Crimper

- High quality construction for durability and long life.
  - Consistent dependable performance.
  - Smooth simple operation.
  - Painted, plated or coated for maximum corrosion resistance
  - Less expensive than competitive models.
  - All parts that contact the container cap are **electrodeless nickel-plated** or **stainless steel**. The benefits include high corrosion resistance and good lubricity have been proven through extensive laboratory and field testing.
  - Textured handle surfaces provide a better grip than smooth painted models. The powder coating used to paint the handles is durable and scratch resistant. The coating process is environmentally clean—no toxic fumes or hazardous wastes are generated. The handle clip, pins and screws are made of corrosion resistant stainless steel.
- Sizes** : 20mm, 13,, 11mm, 8mm for Standard Caps and Plastic Flipoff Caps.


## Vial Caps

- 13mm, 20mm  
Color Coded
- CentreTearOut  
Option Side TearOut  
( to order )
- Plastic FlipOff



## Serum Stoppers

- 13mm, 20mm
- 2-Legged Freeze
- Chlorobutyl Rubber
- Butyl Rubber
- Silicone Rubber
- Teflon Coated



**Vial Kit** : Slab of Vials + Crimper + Caps + Stoppers SPECIAL Discount - **15% off List Price**  
**Or if Order value >\$500**

## Vial Accessories

- Syringe Filters
- Filter Funnels, Flasks  
( & Vacuum Pump )
- Membrane Filters
- Stoppers
- Caps
- Syringe Dispenser
- Hand Crimper  
also
- Bench Crimper
- Pneumatic Crimper



<b>Crimp Vial KITS</b>	More PDF Details <a href="#">PDF</a> <b>Prices</b> Ask for PDF Price List	<a href="#">Lab Dispensers</a>	Filters : Syringe Filters Membrane Filters, VacuumFiltration Equipment
------------------------	---	--------------------------------	--

**for Convenience**  
**Vial Kits**  
 - Vials 1 x slab ( *see Slab Qty* )  
 - Clear or Amber  
 - Volume 10, 13, or 20ml  
 - Vial Top 20mm  
 Chlorobutyl Stoppers RB *matching Qty*  
 Crimp Caps - Color(X) - Standard, Flip-Off Kimble Or WestPharmaceutical (embossed) Flip Off *matching Qty*  
 Hand Crimper - *each Size to match Vials and Caps*

**Discount  
15%  
off  
List Price**

To Order Kit  
 CT Vial Kit Clear or Amber      Vol(ml)    20mm Top Stopper    Crimp Cap / Color / Standard / FlipOff Type/Kimble/WP  
 CT VK-                                    C / A                                    10 / 13 / 20V    20T    -    CB                                    C                                    ( X )                                    S                                    KFO / WPFO  
 (X) specify Color

Matching . . .  
 Vial Slab  
 Stoppers - same Qty  
 Crimp Cap - same Qty  
 Hand Crimper

VIALS in Slabs Shrink-wrapped Clean - NOT Sterile	VIALS Cat.#	Qty/Slab	Kg/Slab
Clear	CTSC-2V-13T-15D32H	371	1.5
	CTSC-4V-13T-16D38H	368	2.0
	CTSC-7V-20T-22D38H	150	1.3
	CTSC-10V-20T-22D45H *	198*	1.9
	CTSC-10V-13T-22D45H	198	1.8
	CTSC-13V-20T-22D49H *	155*	1.9
	CTSC-13V-13T-22D50H	165	1.9
	CTSC-20V-20T-26D66H	145	2.2
	CTSC-23V-20T-29D60H	85	1.4
	CTSCW-50V-20T-43D73H	112	6.3
CTSCW-100V-20T-52D95H	66	6.7	
Amber	CTSA-2V-13T-15D32H	371	1.7
	CTSA-7V-20T-22D38H	241	2.2
	CTSA-13V-20T-24D50HTP *	155*	2.1
	CTSA-23V-20T-29D60HTP	113	2.2
	CTSAW-30V-20T-36D63H	92	7
	CTSAW-50V-20T-43D73H	112	7.4
	CTSAW-100V-20T-52D94H	66	-
Items in Blue - most popular Kits - individual items available see <b>Price List - ASK!</b> * \$1.08 per Set 3 parts			



Clear Vials



50/20      23/20      20/20      13/20  
4/13      2/13      13/10      10/20      10/13      7/20

*Nominal Volume (ml) / actual Liquid Volume / Neck Size(mm) - Stoppers shown in place No Alum Caps shown*

*Vials supplied Clean - Shrink-wrapped in slabs (NOT Sterile) and Empty*

Amber Vials



100/20

23/20

13/20

Stoppers 20mm ( and 13mm )



CT-20BNRB-S

CT-20BNRB-2L

Round Bottom

2-Leg Freeze Dry

Standard

Custom

