



Cat.# 22358

1

Electronic Rechargeable Crimpers and Decappers

- Easy to use; comfortable grip.
- Hundreds of operations from one charge.
- Adjustable crimping force.
- Large battery allows for faster charge cycle.
- All kits include a universal plug set for the battery charger.
- For 11 mm and 20 mm caps.
- One-year warranty.

Instruction Sheets



Select a Product

Catalog #	Product Name	Size	Units	
22358	Electronic Crimper	11 mm	kit	Select
22359	Electronic Crimper	20 mm	kit	Select
22360	Electronic Decapper	11 mm	kit	Select
22361	Electronic Decapper	20 mm	kit	Select

Showing 1 to 4 of 4

Description

- Easy to use; comfortable grip.
- Hundreds of operations from one charge.
- Adjustable crimping force.
- Large battery allows for faster charge cycle.
- All kits include a universal plug set for the battery charger.
- For 11 mm and 20 mm caps.
- One-year warranty.

These electronic crimpers and decappers have a rechargeable 6.4-volt battery installed. The large battery allows for quick, powerful operation and a faster charge cycle, minimizing downtime. Charging the battery takes thirty minutes to two hours. The universal plug set includes U.S., European, UK, and Australian plug clips.