

Disposable Syringes & Kits



New Era supplies an assortment of disposable plastic syringes between 1mL to 140mL from Monoject, NormJect, Terumo & Excel. You can receive discounts in box quantities.

| Description (mL=CC) | Part # | Price Each | Box Part # | Box Price |
|----------------------------|-------------|------------|---------------|-----------------------|
| 140 mL Luer Lock Non-Inert | SYR-MJ140LL | \$6.00 | SYR-MJ140LLBX | \$119.00 (20 per Box) |
| 60 mL Luer Lock Non-Inert | SYR-TER60LL | \$1.25 | SYR-TER60LLBX | \$18.75 (25 per Box) |
| 50 mL Luer Lock Inert | SYR-NJ50LL | \$1.75 | SYR-NJ50LLBX | \$37.50 (30 per Box) |
| 30 mL Luer Lock Inert | SYR-NJ30LL | \$1.25 | SYR-NJ30LLBX | \$40.00 (50 per Box) |
| 20 mL Luer Lock Inert | SYR-NJ20LL | \$0.75 | SYR-NJ20LLBX | \$55.00 (100 per Box) |
| 10 mL Luer Lock Inert | SYR-NJ10LL | \$0.50 | SYR-NJ10LLBX | \$40.00 (100 per Box) |
| 5 mL Luer Lock Inert | SYR-NJ5LL | \$0.40 | SYR-NJ5LLBX | \$30.00 (100 per Box) |
| 3 mL Luer Lock Inert | SYR-NJ3LL | \$0.35 | SYR-NJ3LLBX | \$25.00 (100 per Box) |
| 1 mL Luer-Lock Non-Inert | SYR-EX1LL | \$0.40 | SYR-EX1LLBX | \$30.00 (100 per Box) |
| 1 mL Tuberculin Inert | SYR-TB1 | \$0.25 | SYR-EX1LLBX | \$20.00 (100 per Box) |
| Large Syringe kit | P-SYRKIT-LG | \$15 | N/A | N/A |
| Small Syringe kit | P-SYRKIT-SM | \$11.50 | N/A | N/A |
| Dual Pump Plumbing Kit | P-DKIT | \$20.00 | N/A | N/A |



Large Syringe Kit
Part ID: P-SYRKIT-LG

- (2) 20mL Luer Lock Inert
- (2) 30mL Luer Lock Inert
- (2) 60mL Terumo Luer Lock Non-Inert
- (6) Female Luer Barbs
- (3) Male Luer Barbs
- (20 ft) Tubing .125" ID; .188" OD, PVC



Small Syringe Kit
Part ID: P-SYRKIT-SM

- (2) 1mL Luer Lock Inert
- (2) 3mL Luer Lock Inert
- (2) 5mL Luer Lock Inert
- (2) 10mL Luer Lock Inert
- (6) Female Luer Barbs
- (3) Male Luer Barbs
- (20 ft) Tubing .125" ID; .188" OD, PVC



Dual Pump Plumbing Kit
Part ID: P-DKIT

- For use with our continuous infusion system or dual infusion system (Dual-NE-1000). Provides a dual check valve for two syringes, withdraw will pull fluid from a reservoir and infusion will output into a separate line
- The kit is pre-assembled, all you need to do is plug in the syringes
- Each dual check valve uses 6" of 1/8" tubing to connect to a Y connector, one for reservoir and one for output; each Y-connector has 5' of tubing to make connections to your application. Two 60cc Terumo syringes are included

Highly corrosive chemicals are not compatible with this kit. Contact us for a solution that is compatible.



SyringePump.com
Clever Pumps, Priced Right!

138 Toledo St.
Farmingdale, NY 11735
Phone: (631) 249-1392
Fax: (707) 248-2089
info@syringepump.com