

E-Fill

Electronic Pipette Filling Device

Easy
Efficient
Electronic
Economical

Features

- Two LEDs to show charging status and low battery.
- Continuous adjustment of pump speed during operations.
- An universal Silicon collet is provided to accept all standard pipettes up to 100ml.
- Internal knurling on the collet holds the pipette firmly.
- Comfortable grip.
- Line membrane filter of 0.2 micron prevents Aerosol contamination and liquid from entering the equipment by accident.
- Ni-MH 3.6 rechargeable battery allows eight hours of continuous working.

Microlit



Ordering Information and Specifications

Model No.	Description
MPF - GY - 220V	Grey body/ yellow collet and plungers
MPF - GY - 110V	Grey body/ yellow collet and plungers
MPF - GR - 220V	Grey body/ red collet and plungers
MPF - GR - 110V	Grey body/ red collet and plungers
MPF - GV - 220V	Grey body/ Violet collet and plungers
MPF - GV - 110V	Grey body/ Violet collet and plungers