

Electronic Micropipettes



FEATURES

- Based on Stepper Motor Technology.
- User friendly software.
- Recharging during operation.
- Optimized ergonomics and working comfort.
- Built-in tip ejector.
- Six models cover the complete Pipetting range from 0.2 μl to 10000 μl .



Mode (Easy selection of modes from single key)

- Stepper mode
- Dilution mode
- Standard mode
- Calibration mode

Variable Working Speed can be selected on a scale of 1-5

Easy volume selection from 'UP' & 'DOWN' keys.

Large Fire Button

Single 'FIRE BUTTON' for filling and dispensing.

Large Liquid Crystal Display

Two line Alphanumeric LCD.

Calibration

User friendly calibration mode facilitates easy calibration by the user.

Ordering Information and Specifications

Model	Vol. Range	Increments	Tip Style
ME - 10	0.2-10 μl	0.05 μl	10 μl
ME - 120	5-120 μl	0.5 μl	200 μl
ME - 300	20-300 μl	1.0 μl	300 μl
ME - 1000	50-1000 μl	5.0 μl	1000 μl
ME - 5000	100-5000 μl	10.0 μl	5000 μl
ME - 10000	500-10000 μl	50.0 μl	10000 μl