

# Fluode Gas Sampling Bag - *The NEW Standard in Gas Sampling*

Fluode Gas Sampling Bags Fluode Gas Sampling Bag with stopcock side-opening valve with silicone septum port 1/4" 6.35m FLD21\_0.5L

## ALL PRODUCTS

Tedlar® PVF Gas Sampling Bags

Aluminium-foil (Multi-Layer) Gas Sampling Bags

Teflon® FEP Gas Sampling Bags

Gas Detector/ Monitor—GASTEC/ HONEYWELL/ Riken Keiki/ QRAE

DEVEX (Multi-Layer) Gas Sampling Bags

Kynar PVDF Gas Sampling Bags

New DEVEX (Multi-Layer) Gas Sampling Bags

EVOH Gas Sampling Bags

Polyester Gas Sampling Bags(Odor Bag)

Fluode Gas Sampling Bags

Valves & Fittings for Gas Sampling Bags

## Fluode Gas Sampling Bag with stopcock side-opening valve with silicone septum port 1/4" 6.35m FLD21\_0.5L



**Large Image** : Fluode Gas Sampling Bag with stopcock side-opening valve with silicone septum port 1/4" 6.35m FLD21\_0.5L

### Product Details:

Place of Origin: **Here !**

Brand Name: **HD**

Model Number: **FLD21\_0.5L**

### Payment & Shipping Terms:

Minimum Order Quantity: **10 Pieces**

Price: **AUD \$12.00 /piece**

Packaging Details: **customized**

Payment Terms: **T/T**

Supply Ability: **100000 Pieces per Month**

[Contact Now](#)

[Chat Now](#)

## Fluode Gas Sampling Bag with stopcock side-opening valve with silicone septum FLD21\_1L (air sample bag) port 1/4" 6.35m



Fluode polyfluo-film is developed by our company Dalian Hede technologies for gas sampling bags. The materials of making Fluode polyfluo film are supplied by a renowned fluoroproducts manufacturer. Fluode is a kind of unique polyfluoro-film with excellent chemical stability and superb physical property. The aggregate form of fluorine atoms in Fluode polymer is similar to that of polyvinyl fluoride such as Tedlar, elasticity modulus and stiffness is also close to PVF, but elongation at break increase, which shows up impact resistance performance of Fluode film, that is Fluode film is more tenacious and not susceptible to fracture.. However Fluode has greater f-atom content in CC bond than other fluoro polymers, so the chemical property and physical property of Fluode are improved significantly. Fluode possesses excellent inertness, low permeability low absorptivity and high strength. Fluode gas sampling bags are suitable for low concentration analysis and VOCs assessment. It is excellent substitute of Tedlar bags. The heat-resistance temperature of Fluode fluorine polymer is 150°C ~ 170°C, easy to heat binding sealed, which is very conducive to being produced gas sampling bags, and Fluode gas sampling bag will be developed into new gas sampling bags industry standards.



## Features:

- superb strength and low permeability and low absorptivity
- excellent chemical inertness and stability

### •Standard Fluode Gas Sampling Bag

**Sizes:0.5L,1L,2L,3L,4L,5L,8L,10L,15L,20L,30L,and 40L 6000L or customized size.**

Select from a variety of sizes to suit your needs.

Specification	
Product Type	Sampling Bag
Material	Fluode
Volume(Liters)	0.5
Model	FLD21_0.5L
Brand	HD

## Applicable scope:

Fluode gas sampling bags can be used in high precision of PPM and PPB level sample analysis,can sample and save all kinds of highly corrosive and chemical activegaseous,liquid samples.

**E-mail:sales@chromtech.net.au**

**Whats UP!App:+61 3 9762 2034**

**Website: <https://www.chromtech.net.au/>**

## Packaging & Shipping:

- Reply your inquiry within 24 hours
- Experienced staffs answer all your questions in fluent English>'Strine !
- Exclusive and unique solution can be provided to our customer by our well-trained and professional engineers and staffs . . . *and just maybe beyond the usual "Norm" !*



## Other Products



Fluode gas sampling bag with side-opening PTFE valve (FLD41\_10L) air sample bags



Fluode gas sampling bag with side-opening PTFE valve+ PTFE fitting (FLD4Z\_4L) air sample bags



Fluode gas sampling bag with PTFE valve and septum valve (FLU31\_1L) air sample bags



Fluode air/gas sampling bag with PTFE straight valve+PTFE fitting silicone septum (FLD3Z 1L) syringe



Fluode Gas Sampling Bag with stopcock combination valve with silicone septum (FLD11 0.5L) air sample

### Teflon® FEP Gas

#### Sampling Bags

Teflon® FEP gas sampling bag with PTFE valve & septum port syringe

### New DEVEX

#### (Multi-Layer) Gas

New DEVEX multi-layer foil air/gas sampling bags with stopcock straight valve and septum valve syringe

### Tedlar® PVF Gas

#### Sampling Bags

Tedlar® PVF Gas Sampling Bag with dual-valve with silicone septum port TDLC32\_200L Tedlar PVF sample VOCs detection

