


Permeation Tubes

VICI Metronics Inc.



G-Cal Devices



G-Calibrator



Dynacalibrators



Gas Standards

2nd Source TO-14A/TO-15 Gas Calibration Standards

- Standards from TWO manufacturers provide second source on one order.
- 12 month stability in transportable cylinders.
- Drop shipped for fast delivery and maximum shelf life.



A. Spectra (Linde)
104L Cylinders

B. Scotty (Air Liquide)
110L Cylinders
(Pi-marked Cylinders
for EU Regulations)

For regulators,



RESTEK

DCG Partnership Cylinders:
Size: 7.6 x 24 cm
CGA-170/110 connection.
US DOT Specs: DOT-4B-240ET



Mini-Regulator



22032

Scotty® (Air Liquide) 14
Contents: 14 liters
Pressure: 240 psig (17 bar)
Outlet Fitting: CGA 160
Weight: 1.5 lbs/0.7 kg
Dimensions: 3" diameter x
11" height (7.6 x 28 cm)
DOT Specifications: 4B240



Scotty® (Air Liquide) 48
Contents: 48 liters
Pressure: 300 psig (21 bar)
Outlet Fitting: CGA 165
Weight: 1.75 lbs/0.8 kg
Dimensions: 4" diameter x
16 1/2" height (10.2 x 41 cm)
DOT Specifications: 39 NRC



Spectra Gas / Linde
Scott Gas / Air Liquide
DCG-Gases