



Home » [Have You Got Cables? – Start Here](#)

Have You Got Cables? – Start Here

If you have got a pump and a PC you will need cables to connect them. Have you got cables?

A lot of folks either expect the cables to come with a new pump or they have come to a place where the cables have been lost – I mean put in a special place for safety :-}

This page will take you to other pages which will help you get cables through purchasing information or DIY cables information.

List of Cables, Application and Part Numbers



CBL-PC-PUMP

This page contains images of all the official pump related cables. It lists when to use each cable (application) and their official part numbers. You can order these cables from your [pump distributor](#).

CBL-PC-PUMP is the cable used to connect a pump to your computer.

CBL-NET is the cable to connect a daisy chain from the first pump to the second pump, pump n to pump (n+1).

Where to buy cables?


The quick answer is from your pump distributor. [Find your pump distributor](#).

Hey tell 'em SyringePumpPro sent you.

Don't know what to order?

The part numbers are on [this diagram](#).

The **CBL-PC** pictured here should definitely be on your shopping list along with a USB-RS232 adapter. [Here is the manual for that cable](#).



Thinking of simply buying some generic RS-232 cables? _____

Don't. Take care I regularly see folks struggling with bought cables and then purchasing the pump manufacturer cables in the end. I have photos of a the handy-work of one frustrated general serial cable purchaser – it reveals hacksaw cutting on a cable end – ugly!

DIY cables? _____



tools

Bought cables are inexpensive and will save you time, simply because they will work properly and leave you in no doubt that the cables are the correct ones.

However a lot of people like to build their own. [Here's how...](#)

What are all the different cables/ what are they used for? _____

There's a [table showing each cable and describing it's use](#) here.

Do you need a power supply for your pump? _____

[Here's how.](#)

[Back to Software Support](#) | [Back to FAQs](#)

About the Author: SyringePumpPro



As the author of SyringePumpPro products I have been involved with laboratory pumps for about 10 years now. My career spans electronics, avionics, programming, teaching, research and development laboratory experience, and even television.

SyringePumpPro

Take control of your pumps!

Search

RECENT POSTS

Pump Programming Triggers

Meet the new look New Era Pump Family

New Era Pump Systems Inc

Parts Of A Syringe Pump

Front Panel Controls

MENU

Home

Tour

Compatible Pumps

Pump Manufacturers

Contact Us

Who's Using

Testimonials

About

Blog

Papers Citing SyringePumpPro

Shop

Pricing

Try SyringePumpPro

Find a Distributor

Licensing

Pump Programming Service

Used Pumps

Login

Pumps

Programming

Syringes and Tubing

Manuals

Rates and Specifications

Drawings

Connectivity

Cables

Applications

TAGS

Aladdin

Braintree Scientific

cables

Chemyx

Cole-Parmer

distributor

diy dual

Electro-*

Harvard

KDScientific

licence

microscope

NE-300

NE-500

NE-501L

NE-510

NE-1000

NE-1010

See full tags →

© Copyright 2019 SyringePumpPro | All Rights Reserved | [Terms of Use](#) | [Privacy Policy](#)

