

CE

Microlit
ISO 9001 : 2008 CERTIFIED

Li'lpet

FIXED VOLUME
MICROPIPETTE

LOW PRICE

HIGH PRECISION

HIGH ACCRACY

2 STEP PLUNGER

UNIQUE TIP CONE

FULLY AUTOCLAVABLE

Technical Specifications :

Model	Capacity	±%Accuracy	± % CV
MM - 5	5 µl	1.5	1.0
MM - 10	10 µl	1.0	1.0
MM - 20	20 µl	0.5	0.5
MM - 25	25 µl	0.5	0.5
MM - 50	50 µl	0.4	0.3
MM - 100	100 µl	0.3	0.3
MM - 200	200 µl	0.5	0.5
MM - 250	250 µl	0.5	0.5
MM - 500	500 µl	0.5	0.5
MM - 1000	1000 µl	0.5	0.5

Features :

- ✓ **OPTIMUM SIZE**
Microlit Li'lpet is 130 mm in length which is optimum size to facilitate the user in handling the Microlit Li'lpet like a standard pipette and not like a syringe
- ✓ **COLOR-CODED**
Color-coded for various capacities for easy identification in laboratory environment
- ✓ **UNIQUE TIP CONE**
Unique tip cone design to accept both micro tip (200µl) and ultra micro tip in pipettes upto 20 µl.
- ✓ **TWO-STAGE PLUNGER**
Two-stage plunger movement ensures high accuracy and precision in pipetting technique. This also facilitates last drop dispensing.
- ✓ **FULLY AUTOCLAVABLE**
Microlit Li'lpet is fully autoclavable at 121°C & 15 psi for sensitive application.

Microlit

629, Pakramau Kursi Road ,
Lucknow 226026
Tel : 0522-2890767, 2890768, 2890761
Fax : 0522-2890760
Email : ajain@microlit.com
Website : www.microlit.com