

## Li'lpet

# Microlit

- Microlit Li'lpet is 130mm in length which is optimum size to facilitate the user in handling the Microlit Li'lpet like a standard pipette and not like a syringe.
- Unique Tip-Cone design to accept both micro tip (200 $\mu$ l) and ultra micro tip ( 25  $\mu$ l ).
- Two stage plunger movement ensures high accuracy and precision in pipetting technique. This also facilitates last drop dispensing.
- Microlit Li'lpet is FULLY AUTOCLAVABLE at 121°C and 15 psi for sensitive applications.
- Colour coded for various capacities for easy identification in laboratory environment.

### Ordering Information and Specifications

Model	Capacity	Accuracy	CV
MM-5	5 $\mu$ l	1.5%	1.0%
MM-10	10 $\mu$ l	1.0%	1.0%
MM-20	20 $\mu$ l	0.5%	0.5%
MM-25	25 $\mu$ l	0.5%	0.5%
MM-50	50 $\mu$ l	0.4%	0.3%
MM-100	100 $\mu$ l	0.3%	0.3%

Model	Capacity	Accuracy	CV
MM-200	200 $\mu$ l	0.5%	0.5%
MM-250	250 $\mu$ l	0.5%	0.5%
MM-500	500 $\mu$ l	0.5%	0.5%
MM-1000	1000 $\mu$ l	0.5%	0.5%

