

A Autoclaves

Model: STE-Ten, Mini-Fast Autoclave



STE-Ten

Voltage V/Hz	Power W	Capacity L	Chamber size Ømm	Packing Size mm	Weight Kg
220/50	1100	1	90x160	570x470x500	30(net:28)

It adopts B-wclass (prED-13060), with 3-times pre-vacuum preceding vacuum drying. Residual humidity is <0.2, vacuum measurement reaches -0.7bar. It sterilizes and sanitizes miscellaneous packed or unpacked medical instruments, whether solid or A-type hollow, whether with multi-apertures or inset pipes. It efficiently rid air from every hollow area of any hollow instrument, such as hand-piece, making sure that water steam reaches every corner. It is installed with BOWIC & DICK which measures the penetration of water steam. Mini printer can be attached to record the process of sterilization.

Model: STE-23L, Class B 23 Liter Autoclave



STE-23L

Voltage V/Hz	Power W	Capacity L	Chamber size Ømm	Packing Size mm	Weight Kg
220/50	2300	23	247x450	760x550x495	64(net:50)

It adopts B-wclass (prED-13060), with 3-times pre-vacuum preceding vacuum drying. Residual humidity is <0.2, vacuum measurement reaches -0.7bar. It sterilizes & sanitizes miscellaneous packed or unpacked medical instruments, whether solid or A-type hollow, whether with multi-apertures or inset pipes. It efficiently rid air from every hollow area of any hollow instrument, such as hand-piece, making sure that water steam reaches every corner. It is installed with BOWIC & DICK which measures the penetration of water steam. Mini printer can be attached to record the process of sterilization.

Models: UTFB/UTFBS/UTKBS, Vertical Pressure Large Autoclave, 30,50,75,100,150,200 Liter

- Computer control auto recycle sterilization
- LCD screen indicate working state
- Temp. scope 50-126 deg C
- Time scope 0-99h
- Overpressure auto-discharging 0.145-0.165MPa
- Automatic shut down with beep reminding after sterilization



UTFB

Stainless steel sterilizing chamber



UTFBS

Full stainless steel



UTKBS

- Full stainless steel
- Hand wheel type of quick-open door structure. Pressure safety lock

Specifications:

MODEL	UTxxx*-30LV	UTxxx*-50LV	UTKB-75LV	UTKBS-100LV	UTKBS-150LV	UTKBS-200LV
VOLUME	30L	50L	75L	100L	150L	200L
POWER SOURCE	220V	220V	220/380V	380V	380V	380V
CONSUMED POWER	3kw	3.5kw	3.5/4.5kw	4.5kw	6.0kw	6.0kw
DIMENSION OF STERILIZING CHAMBER (mm)	Ø350x330	Ø350x520	Ø388x650	Ø500x550	Ø500x760	Ø500x1000
WEIGHT	60kg	80kg	90kg	100kg	120kg	140kg
DIMENSION (mm)	580x580x1000	600x580x1250	600x600x1400	630x630x1300	630x630x1450	630x630x1690

*xxx=FBS,FB,KB,KBS