

W Water Baths, Oil Bath

Orbital Large Shaking Water Bath

Speed	10~250rpm
Orbit diameter	25mm
Speed readout	Tachometer
Temperature control	PID
Temperature readout	Digital
Working temperature range	Ambient+ 5°C ~75°C
Stability	±0.2°C
Power watts	1200W
Optional accessory	Gable Cover
Platform	450X290
Inner tank size (mm)	560X405X200
Size (mm)	740X530X450
Weight (kg) (N.W. / G.W.)	50.0/53.0

Power Supply: AC 110V 60Hz (220V 50/60Hz available)

Number Of Clamps	35	24	15	12	6	2	1
Capacity Of Flask	50ml	125ml	250ml	500ml	1000ml	2000ml	3000ml

- Ideal for a wide range of laboratory applications requiring controlled temp. with orbital shaking motion.
- Suitable for long time operation up to weeks and months.
- Separate switches for speed ON/OFF & heater ON/OFF.
- Individual readout for RPM & digital readout for temp.
- All stainless steel inner tank, platform & clamps.
- Maintenance-free induction brushless motor.
- Differential clamps available.
- Microprocessor based on PID temperature control.
- Speed and temperature adjustable.
- Quality certification: ISO 9001



Model: BS-1836, Boiling Sterilizer



Features:

- All hydraulically pressed 18.8 stainless steel, can withstand electrolysis without corrosion.
- Long life immersion stainless nickel pipe heater & corrosion proof.
- Precision mechanical timer with bell switch control.
- Waterless overheating prevention system with auto power cut-off.
- Build-in low level cut-off device.

Specifications:


- Timer: 0-60 min. mechanical timer with bell
- Heater: 1500 watts
- Thermostat: Temp. over to 150°C auto cut-off power
- Body material: 18.8 stainless steel
- Tray Dimen.: W350xD165xH30(mm)
- Overall Dimension: W525xD190xH170(mm)
- Net/Gross Weight: 3.5/4.0Kgs
- Power Supply: 220V 50-60Hz -7Amp

Model: B-501, Heated Circulating Oil Bath

3 Models with Choice of Temperature Range and Bath Capacity ● Over Heating Safety Device is Equipped ● Wide Range of Temperature Control ● Digital Temperature Setting & Readout ● External Circulating Pump is Available on Special Order

Specifications:

Model	B-501	B-502	B-503
Temp. Range	Room Temp.+5°C to 200 °C	Room Temp.+5°C to 200 °C	Room Temp.+5°C to 300 °C
Stability	±0.1C	±0.1C	±0.1C
Setting/readout	Digital	Digital	Digital
Heater	1.2KW	1.2KW	1.2KW
Stirring	By Built-in Circulation Pump		
Inside Dimension	W19.5xD23.5xH15cm	W39.5xD29.5xH15cm	W30xD30xH20cm
Outside Dimension	W54xD34.5xH35.5cm	W73xD39xH33.5cm	W66xD40xH39.5cm
Volume	6.8L	17.4L	18L
Power Source	110/220V50/60HZ 13A/6.5A	110/220V50/60HZ 13A/6.5A	110/220V50/60HZ 13A/6.5A



Optional With External Circulation Pump