

C Centrifuges-Medium Size

Model: CN-3101, Economy Centrifuge, up to 5000rpm, Max. Capacity 400ml



Features:

LCD speed display • Regressive 30 min. timer or continuous operation • Stable and stepless speed change • Stainless steel inside bowl • Breakdown indication and fuse safety device • Auto current-restraining function • Overheat and overload safety protection • Hidden auto-balance device • Safety door lock and cover open indication • Auto accelerating & decelerating functions. • Auto-cooling air circulation • See-through window • Optional rotor selections for various tube size (See in page 19 in the full product catalog) • Low noise.

Speed(rpm)	Max. RCF(x g)	Max. Capacity(ml)	Timer(min:sec.)	Weight(kg)	W×D×H Size(cm)
500~5,000	3,790	400	0~30 or "ON"	20.5	36×48×34

Model: CN-3301, Brushless General Purpose Centrifuges, up to 6500rpm, Max. Capacity 400ml



Features:

LCD speed display • Regressive 30 min. timer or continuous operation • Quiet brushless motor • Stable and stepless speed change • Stainless steel inside bowl • Breakdown indication & fuse safety device • Auto current-restraining function • Overheat and overload safety protection • Hidden auto-balance device • Safety door lock and cover open indication • Auto accelerating & decelerating functions • Auto-cooling air circulation • See-through window • Optional rotor selections for various tube size • Low noise.

Speed(rpm)	Max. RCF(x g)	Max. Capacity(ml)	Timer(min:sec.)	Weight(kg)	W×D×H Size(cm)
500~6,500	4,500	400	0~30 or "ON"	18	36×48×34

Model: CN-350, Medium Size Centrifuge, up to 6500rpm, Max. Capacity 400ml



Features:

LCD, speed and time display with RCF setting • Quiet brushless motor • Microprocessor control • 10 memory programs for preset parameters • Auto power off when open • Auto cover lock with motor running • Overheat and overload safety protection • Auto power off for rotor imbalance • Auto buzzer when finish • Hidden auto-balance device • Power pilot light and fuse safety device • Auto accelerating & decelerating functions • See-through window • Optional rotor selections for various tube size • Low noise.

Speed(rpm)	Max. RCF(x g)	Max. Capacity(ml)	Timer(min:sec.)	Weight(kg)	W×D×H Size(cm)
500~6,500	4,500	400	99:59	20.5	36×48×34

Model: CN-350RD, Medium Size Refrigerated Centrifuge, up to 5000rpm, Max. Capacity 480ml



Features:

LCD, temperature, speed and RCF setting also display • Temp. range : 0°C~40°C • 10 memory programs for preset parameters • Overheat and overload safety protection • Stainless steel inside bowl • Microprocessor control • Auto power off for rotor imbalance • Auto buzzer when finish • High function & high effective control system • Breakdown indication & fuse safety device • Auto accelerating & decelerating functions • Auto-cooling air circulation • Low noise.

Speed(rpm)	Max. RCF(x g)	Max. Capacity(ml)	Timer(min:sec.)	Weight(kg)	W×D×H Size(cm)
500~5,000	4,640	480	0~30 or "ON"	40	45×56×38