

# C Centrifuges-Large Size

## Model: CN-5101, Economy Centrifuge, Up to 5000rpm, Max Capacity 1000ml



### Features:

LCD speed display • Regressive 30 min. timer or continuous operation • Stable and stepless speed change • Stainless steel inside bowl • Breakdown indication and fuse safety device • Auto current-restraining function • Overheat & overload safety protection • Hidden auto-balance device • Auto cover lock with cover open indication • Auto accelerating & decelerating functions • Auto-cooling air circulation • See-through window. Optional rotor selections for various tube size • Low noise.

Speed	Max.RCF	Max.Capacity	Timer (min)	Weight	W×D×H Size
500~5,000rpm	4,640x g	480ml	0 30 or "ON"	38kg	45×56×38cm

## Model: CN-5301, Economy Centrifuge, Up to 6500rpm, Max Capacity 1000ml



### Features:

LCD speed display • Regressive 30 min. timer or continuous operation • Stable and stepless speed change • Quiet brushless motor • Stainless steel inside bowl • Breakdown indication and fuse safety device • Auto current-restraining function • Overheat & overload safety protection • Hidden auto-balance device • Auto cover lock with cover open indication • Auto accelerating & decelerating functions • Auto-cooling air circulation • See-through window • Optional rotor selections for various tube size • Low noise.

Speed	Max.RCF	Max.Capacity	Timer (min)	Weight	W×D×H Size
500~6,500rpm	4,900x g	1000ml	0 30 or "ON"	38kg	45×56×38cm

## Model: CN-650, Large Size Centrifuge, Up to 6500rpm, Max Capacity 1000ml



### Features:

LCD, speed & time display with RCF setting • Quiet brushless motor • Microprocessor control • 10 memory programs for preset parameters • Auto power off when open • Auto cover lock with motor running • Overheat & overload safety protection • Auto power off for rotor imbalance • Auto buzzer when finish • Hidden auto-balance device • Power pilot light and fuse safety device • Auto accelerating and decelerating functions • See-through window • Optional rotor selections for various tube size • Low noise.

Speed	Max.RCF	Max.Capacity	Timer	Weight	W×D×H Size
500~6,500rpm	4,900x g	1000ml	99:59min:sec.	40kg	45×56×38cm

## Model: CN-650R, Large Size Refrigerated Centrifuge, Up to 6500rpm, Max Capacity 1000ml



### Features:

LCD, temperature, speed and RCF setting and display • Temperature range: 4°C~40°C • 10 memory programs for preset parameters • Overheat & overload safety protection • Stainless steel inside bowl • Microprocessor control • Auto power off for rotor imbalance • Auto buzzer when finish • Auto current-restraining function • Breakdown indication & fuse safety device • Auto accelerating and decelerating functions • Auto-cooling air circulation • Low noise.

Speed	Max.RCF	Max.Capacity	Timer	Weight	W×D×H Size
500~6,500rpm	4,900x g	1000ml	99:59min:sec.	112kg	52×58×90cm

## Model: CN-650RD, Large Size Refrigerated Centrifuge, Up to 6500rpm, Max Capacity 1000ml





























### Features:

LCD, temperature, speed and RCF setting and display • Temperature range: 4°C~40°C • 10 memory programs for preset parameters • Overheat & overload safety protection • Stainless steel inside bowl • Microprocessor control • Auto power off for rotor imbalance. • Auto buzzer when finish • Auto current-restraining function • Breakdown indication & fuse safety device • Auto accelerating & decelerating functions • Auto-cooling air circulation • Low noise.

Speed	Max.RCF	Max.Capacity	Timer	Weight	W×D×H Size
500~6,500rpm	4,900x g	1000ml	99:59min:sec.	75kg	80×56×38cm

# Centrifuges-Large Size C

**Models: CN-650R, CN-650RD, CN-650, CN-5101, CN-5301, Optional Rotors & Adapters For Centrifuges**

Applicable Model	Rotor		Max. Speed (rpm)	Max. RCF (xg)	Capacity (ml x pcs)	Tube Size (φ x l)mm		Cat. No.	
						Max.	Min.		
CN-650R CN-650RD CN-650 CN-5301	RS-250		3,500	2,250	1.5x96	10.7x100	8x41	2-6000-21	
					6x112	12x110	8.5x60	2-6000-22	
					10x80	15.8x100	12x60	2-6000-23	
					15x40TC	16.5x122	14x60	2-6000-25	
					15x48	17.2x122	14.5x60	2-6000-24	
					50x16TC	29.5x125	27.5x50	2-6000-26	
						250x4	62x130	58x70	2-6000-27
									
			<b>2-6000-21</b>	<b>2-6000-22</b>	<b>2-6000-23</b>	<b>2-6000-25</b>	<b>2-6000-24</b>	<b>2-6000-26</b>	<b>2-6000-27</b>
	RS-1548		3,500	2,430	15x48	17.5x107	14x80	2-6000-12	
RS-2504		3,500	2,480	250x4	65.5x130	61x90	2-6000-13		
RS-9602		2,500	1,060	96micro x4	129x86	126x83	2-6000-14		
CN-650R CN-650RD CN-650 CN-5101 CN-5301	RS-100		3,500	2,090	6x56	13x105	10x50	2-5400-21	
			3,500	2,150	10x48	16x107	13x50	2-5400-22	
			3,500	2,150	15x16TC	16.7x135	13.5x86	2-3400-32	
			3,500	2,150	15x32	17.5x115	14x80	2-5400-23	
			4,000	2,770	50x4	35x135	32x90	2-3400-22	
			5,000	4,640	50x4TC	29.5x137	26.5x90	2-3400-23	
			3,500	2,120	50x8	35x120	32x90	2-5400-24	
			3,500	2,120	50x8TC	29.5x127	26.5x90	2-5400-25	
			4,000	2,930	100x4	45.5x130	42x95	2-3400-21	
									
		<b>2-5400-21</b>	<b>2-5400-22</b>	<b>2-3400-32</b>	<b>2-5400-23</b>	<b>2-3400-23</b>	<b>2-5400-24</b>	<b>2-5400-25</b>	<b>2-3400-21</b>
RS-50		4,500	3,350	15x16	16.7x115	13.5x86	2-3400-32		
		4,000	2,600	50x4	35x117	32x90	2-3400-22		
		5,500	5,270	50x4TC	29.5x120	26.5x90	2-3400-23		
		5,000	4,200	100x4	45x112	42x95	2-3400-21		
									
		<b>2-3400-32</b>	<b>2-3400-22</b>	<b>2-3400-23</b>	<b>2-3400-21</b>				
RA-5008		6,000	4,500	50x8	29.5x118	27.5x90	2-6000-16		
RA-1512H		6,000	4,500	15x12	16.7x130	13.5x86	2-6000-15		