

C Centrifuges-Large Size

Model: CN-5101, Economy Centrifuge, Up to 5000rpm, Max Capacity 1000ml



Features:

LCD speed display • Regressive 30 min. timer or continuous operation • Stable and stepless speed change • Stainless steel inside bowl • Breakdown indication and fuse safety device • Auto current-restraining function • Overheat & overload safety protection • Hidden auto-balance device • Auto cover lock with cover open indication • Auto accelerating & decelerating functions • Auto-cooling air circulation • See-through window. Optional rotor selections for various tube size • Low noise.

Speed	Max.RCF	Max.Capacity	Timer (min)	Weight	W×D×H Size
500~5,000rpm	4,640x g	480ml	0 30 or "ON"	38kg	45×56×38cm

Model: CN-5301, Economy Centrifuge, Up to 6500rpm, Max Capacity 1000ml



Features:

LCD speed display • Regressive 30 min. timer or continuous operation • Stable and stepless speed change • Quiet brushless motor • Stainless steel inside bowl • Breakdown indication and fuse safety device • Auto current-restraining function • Overheat & overload safety protection • Hidden auto-balance device • Auto cover lock with cover open indication • Auto accelerating & decelerating functions • Auto-cooling air circulation • See-through window • Optional rotor selections for various tube size • Low noise.

Speed	Max.RCF	Max.Capacity	Timer (min)	Weight	W×D×H Size
500~6,500rpm	4,900x g	1000ml	0 30 or "ON"	38kg	45×56×38cm

Model: CN-650, Large Size Centrifuge, Up to 6500rpm, Max Capacity 1000ml



Features:

LCD, speed & time display with RCF setting • Quiet brushless motor • Microprocessor control • 10 memory programs for preset parameters • Auto power off when open • Auto cover lock with motor running • Overheat & overload safety protection • Auto power off for rotor imbalance • Auto buzzer when finish • Hidden auto-balance device • Power pilot light and fuse safety device • Auto accelerating and decelerating functions • See-through window • Optional rotor selections for various tube size • Low noise.

Speed	Max.RCF	Max.Capacity	Timer	Weight	W×D×H Size
500~6,500rpm	4,900x g	1000ml	99:59min:sec.	40kg	45×56×38cm

Model: CN-650R, Large Size Refrigerated Centrifuge, Up to 6500rpm, Max Capacity 1000ml



Features:

LCD, temperature, speed and RCF setting and display • Temperature range: 4°C~40°C • 10 memory programs for preset parameters • Overheat & overload safety protection • Stainless steel inside bowl • Microprocessor control • Auto power off for rotor imbalance • Auto buzzer when finish • Auto current-restraining function • Breakdown indication & fuse safety device • Auto accelerating and decelerating functions • Auto-cooling air circulation • Low noise.

Speed	Max.RCF	Max.Capacity	Timer	Weight	W×D×H Size
500~6,500rpm	4,900x g	1000ml	99:59min:sec.	112kg	52×58×90cm

Model: CN-650RD, Large Size Refrigerated Centrifuge, Up to 6500rpm, Max Capacity 1000ml



Features:

LCD, temperature, speed and RCF setting and display • Temperature range: 4°C~40°C • 10 memory programs for preset parameters • Overheat & overload safety protection • Stainless steel inside bowl • Microprocessor control • Auto power off for rotor imbalance. • Auto buzzer when finish • Auto current-restraining function • Breakdown indication & fuse safety device • Auto accelerating & decelerating functions • Auto-cooling air circulation • Low noise.

Speed	Max.RCF	Max.Capacity	Timer	Weight	W×D×H Size
500~6,500rpm	4,900x g	1000ml	99:59min:sec.	75kg	80×56×38cm