



Specialty Gases & Equipment



It is the guiding principle of MESA Specialty Gas and equipment to offer only the highest quality calibration gas products to our customers.

Service...

We recognize that only excellent customer service will allow our growth in this highly competitive calibration gas industry.

Support...

MESA Specialty Gas offers technical support through top experts in the field. Let us know how we can assist you.

Natural Gas and BTU Standards

The following BTU mixes are held in inventory and are available for immediate delivery in small disposable cylinders. For cylinder information, please see cylinder specifications for the [20L and 14L Cylinder](#). The recommended regulator for use with the 20L and 14L is a [Series 300](#) regulator.

STOCK #:	K215585	03143	360816	M-03085	99282	X-20987
COMPONENTS	CONCENTRATION	CONCENTRATION	CONCENTRATION	CONCENTRATION	CONCENTRATION	CONCENTRATION
NITROGEN	1.53%	2.39%	0.23%	0.95%	9.11%	0.23%
CARBON DIOXIDE		1.91%	10.05%	6.91%	2.97%	10.44%
METHANE	BALANCE	BALANCE	BALANCE	BALANCE	BALANCE	BALANCE
ETHANE	2.49%	4.80%	4.19%	6.50%	2.94%	9.95%
PROPANE	2.51%	0.973%	1.87%	4.63%	1.00%	4.01%
ISO-BUTANE	0.50%	0.286%	0.880%	0.919%	0.398%	0.890%
N-BUTANE	2.00%	0.288%	0.301%	0.920%	1.010%	0.837%
NEOPENTANE					0.100%	
ISO-PENTANE	0.301%	0.095%	0.210%	0.461%	0.201%	0.398%
N-PENTANE	0.500%	0.095%	0.110%	0.463%	0.397%	0.349%
HEXANE	0.201%	0.077%	0.372%	0.324%	0.050%	0.298%
HEPTANE		0.006%				

LIST IN 14L (USD)

LIST IN 20L (USD)



Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

RGA Mixes

The following RGA mixes are held in inventory and are available for immediate delivery in small disposable cylinders. For cylinder information, please see cylinder specifications for the [20L and 14L Cylinder](#). The recommended regulator for use with the 20L and 14L is a [Series 300](#) regulator.

STOCK #:	47968AW	32286AW	8994AW
COMPONENTS	<u>CONCENTRATION</u>	<u>CONCENTRATION</u>	<u>CONCENTRATION</u>
HELIUM			
HYDROGEN	BALANCE	3.26%	12.01%
NITROGEN	3.91%	BALANCE	BALANCE
OXYGEN			
ARGON	0.49%	0.17%	0.99%
CARBON MONOXIDE	0.98%		1.01%
CARBON DIOXIDE	2.94%	0.10%	3.02%
METHANE	8.23%	1.99%	5.00%
ETHANE	5.84%	0.51%	3.99%
PROPANE	6.83%	0.94%	5.95%
CYCLOPROPANE	0.05%		0.03%
ISO-BUTANE	4.91%	2.07%	5.14%
N-BUTANE	3.86%	0.54%	3.98%
NEOPENTANE			
ISO-PENTANE	0.97%	1.48%	0.99%
N-PENTANE	0.98%	0.49%	0.98%
HEXANE	0.49%	0.59%	0.10%
ETHYLENE	1.97%	0.51%	1.96%
ACETYLENE	1.94%		1.00%
PROPYLENE	1.95%	2.36%	2.99%
PROPADIENE	1.76%		0.99%
1-BUTENE	1.96%	0.62%	2.00%
CIS-2-BUTENE	1.96%	0.90%	0.99%
TRANS-2-BUTENE	1.95%	1.09%	2.96%
ISOBUTYLENE	1.94%	1.05%	0.99%
1,2-BUTADIENE			
1,3-BUTADIENE	2.91%	0.12%	3.04%
METHYL ACETYLENE			
VINYL ACETYLENE			
ETHYL ACETYLENE			
1-PENTENE		0.59%	0.39%
TRANS-2-PENTENE		0.30%	0.20%
CIS-2-PENTENE		0.20%	0.30%
2-METHYL-2-BUTENE		0.41%	0.19%
2-METHYL-1-BUTENE			
3-METHYL-1-BUTENE		0.14%	

3-METHYL-1,2-BUTADIENE

LIST IN 14L (USD)

LIST IN 20L (USD)

HROM*alytic* +61(0)3 9762 2034
ECH*nology* Pty Ltd

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA