



## Custom Gas Chromatographs

### MENU

[Home](#)

[About Us](#)

[Products](#) >>

[Downloads](#) >>

[Blog](#)

[Distributors](#)

[Contact Us](#)

[Privacy Policy](#)

### FREE CONSULTATION

FREE CONSULTATION TO DISCUSS  
YOUR GC APPLICATION 310 214 5092

[Call Us Now](#)

### SUBMIT INQUIRY



REQUEST A  
QUOTE



### INTERACTIVE PRODUCT GUIDE



## MODEL 310 CANNABIS AND HEMP GC SYSTEM

[Home](#) > [Gas chromatographs](#) > Model 310 Cannabis and Hemp GC system



310-MM

### Model 310 Cannabis and Hemp GC system

#### Compact Multi-Purpose Cannabis and Hemp GC Analyzer

The model 310 cannabis and hemp GC system tests cannabis and hemp for cannabinoids (THC, CBD, CBN, etc.), terpenes (pinene, myrcene, limonene, etc.), as well as residual solvents (butane, ethanol, etc.) that may be present in any concentrates.

It is a fast, simple, and precise way to test your product, while providing a laboratory grade analysis. No prior experience with lab equipment is necessary. One half day (four hours) of free training, good for up to four people, is provided at our Los Angeles area offices, which are conveniently located near Los Angeles International Airport (LAX).

No external gas cylinders are required to operate the model 310 Cannabis and Hemp GC system, as it is equipped with a built-in hydrogen generator and air compressor.

The 310 Cannabis and Hemp GC system features:

- a flame ionization detector (FID)
- a built-in hydrogen generator
- a built-in air compressor
- a 30meter MXT-502 capillary column for cannabinoid and terpene measurement
- a Hayesep D packed column for residual solvent measurement
- two on-column injection ports so columns do not have to be exchanged
- a built-in PeakSimple USB data system (includes PeakSimple software for Windows PC at no extra charge)
- free technical support and software updates for the life of the GC system

(Prices as of January 1, 2024. PRICES ARE SUBJECT TO CHANGE. CONSULT MOST RECENT PRICELIST)

0310-0091	Model 310 Cannabis and Hemp GC system, 115VAC	\$11,926.00
0310-0091-2	Model 310 Cannabis and Hemp GC system, 230VAC	\$11,926.00
0310-0095	Model 310 Cannabis Edibles GC system, 115VAC	\$13,814.00
0310-0095-2	Model 310 Cannabis Edibles GC system, 230VAC	\$13,814.00

### Select Model

--Select Model--



[Send Inquiry](#)

or

[Add To Inquiry](#)

### About SRI Instruments

SRI Instruments is a manufacturer of high quality chromatography instruments and accessories. Our GC systems offer fast and simple cannabis analysis solutions. We take pride in providing our customers with exceptional service and attention, which is reflected in our high level of customer satisfaction. When you call SRI Instruments, you get the answers you seek from our knowledgeable technical staff.



[Home](#) | [About Us](#) | [Blogs](#) | [Contact Us](#)