


maybe ? even **brilliant** ! for Education Systems > kindergarten UP !  
 > a "FullBlown" PC > GC Data System ( PeakSimple )

## Model 30 TCD Mini Gas Chromatograph and PeakSimple Data System



Complete **GC** and **Data System**  
 Super Portable ( book size )  
 Ambient temperature operation  
 No gas cylinders-runs on room air  
 Packed or capillary columns  
 Includes 2 meter packed column  
 TCD detector PeakSimple USB  
 Data System Made in USA  
 -Two year warranty  
 Low cost (USD) **\$3,687.00**





Vacuum pump

The SRI Instruments Model 30 Gas Chromatograph is a simple low-cost system intended for basic GC applications such as measuring Hydrogen in Air or alcohol in water. The extremely small size, USB connection to any Windows computer ( laptop ) and ability to operate on room air using the included vacuum pump makes the Model 30 perfect for field operations. Helium carrier gas may also be used with the optional carrier gas pressure regulator kit.

(2022 pricing, prices subject to change, consult most recent price list.)


### Order:

0030-1003 Model 30 Basic TCD GC \$3,687.00  
 8670-1358 Pressure Regulator w/ Gauge \$513.00



TCD detector

column



On-column injector

email: [sales@srigc.com](mailto:sales@srigc.com) > > ? Check but **Be-A-WARE** ! > your learning curve !