

# NE-300 "Just Infusion"™ Syringe Pump



<b>Models</b>
NE-1000 Single Syringe Pump
NE-4000 Double Syringe Pump
NE-1600/NE-1800 Multi-Syringe Pumps
NE-500 OEM Syringe Pump
NE-300 "Just Infusion"™
NE-1002X Microfluidics Syringe Pump
NE-9000 Dispensing Peristaltic Pump
NE-7000 Adhesive Dispenser
Disruptor, Homogenizer & Rocker
<b>Accessories</b>
Syringes and Plumbing Supplies
Pump Accessories
Syringe Heater
Closed Loop Analog Sensor Interface

## NE-300 "Just Infusion"™ Syringe Pump



**\$275 Breakthrough Pricing!!!**

NE-300 "Just Infusion"™ Syringe Pump  
NE-300

Price: USD \$275

Shipping Charges are added onto order, Domestic orders only.  
International orders please contact us for more information

("Add to cart" button [Google Checkout] is for buying the NE-300 only.)



Australian Distributors  
Importers & Manufacturers  
www.chromtech.net.au

11/12

## [Special Applications](#)

[Continuous Infusion,  
Dual Pump System](#)

[Linear/Gradient  
Pumping Upgrade](#)

[Pumping Program  
Memory Upgrade](#)

[OEM Pump & Other  
Design Services](#)

## [Resources](#)

[Research Papers](#)

[Frequently Asked  
Questions](#)

## Description:

"Just the features that you want, and none that you don't want to pay for."

- Holds small or big syringes up to 60 mL
- Pumps continuously until you stop the pump
- Infusion rate can be changed while pumping
- Remembers previous settings on power up
- Power failure mode allows pump to continue pumping after a reset
- Syringe purge mode
- Displays total volume dispensed in mL or  $\mu$ L units
- Maximum pumping rate of 1500 mL/hr with a 60 mL syringe
- Precise reproducible flow rates
- Selectable infusion rates units: mL/hr,  $\mu$ L/hr, mL/min,  $\mu$ L/min
- Pictured Syringe is for Demonstration purposes only.
- [Disposable, Glass, Stainless Steel Syringes and Plumbing Supplies](#) are sold separately.

"Infuses Only!"

- -Does not withdraw, volume target, have computer interface or have programmability
- -Just the features that you want to pay for
- - [NE-300 Manual](#)

[International Orders \(click here for important details\)](#)

\*\*\* Not for clinical use on humans \*\*\*



Australian Distributors  
Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)

11/12