

Stainless Steel Syringes

Maximum & Minimum Flow Rates

NE-1000*

Syringe (cc)	Inside Diameter (mm)	Maximum Rate (mL/hr)	Minimum Rate (µl/hr)	Maximum Rate (mL/min)
1ml	9.538	218.6	3.005	3.644
3ml	9.538	218.6	3.005	3.644
5ml	12.7	387.6	5.326	6.46
8ml	9.538	218.6	3.005	3.644
20ml	19.13	879.5	12.09	14.65
50ml	28.6	1965	27.01	32.76
100ml	34.93	2932	40.29	48.87
200ml	44.75	4813	66.13	80.21
1ml	9.538	787.5	6.009	13.12
3ml	9.538	787.5	6.009	13.12
5ml	12.7	1396	10.66	23.26
8ml	9.538	787.5	6.009	13.12
20ml	19.13	3167	24.17	52.79
50ml	28.6	7080	54.03	118
100ml	34.93	9999	80.58	176
200ml	44.75	9999	132.3	288.9

NE-1010**

*Flow rates for the NE-1000 also apply to the following pumps: NE-500, NE-501, NE-500L, & NE-501L.

**Flow rates for the NE-1010 also apply to the following pumps: NE-510, NE-511, NE-510L, & NE-511L.