

# Netbook Syringe Pump - \$995



## Features of the Netbook:

- Intel Atom™ Processor
- Windows 7 Starter
- 10.1" Screen
- Wireless Internet Adapter

## Features of the NE-500 Syringe Pump

- Infusion and withdrawal
- Set a single pumping rate and/or dispensing volume
- Holds 1 syringe up to 60 cc
- Infusion rates from 0.73  $\mu\text{L/hr}$  ( 1 cc syringe) to 2100 ml/hr (60 cc syringe)
- Dispensing accuracy of  $\pm 1\%$
- Program up to 41 pumping phases that change pumping rates, set dispensing volumes, insert pauses, control and respond to external signals, sound the buzzer
- Includes all the advanced functions and programming features of the NE-1000 family
- Disposable, Glass, Stainless Steel Syringes and Plumbing Supplies are sold separately.
- Unlimited lifetime technical support
- Two year warranty

## Description:

- User friendly netbook computer control of syringe pump
- Comes with all necessary components for full standalone usage out of the box

## Includes:

- Basic Netbook Computer
- License copy of Syringe Pump Pro software
- NE-500 OEM Syringe Pump
- Small Syringe Kit
- RS-232 To USB Converter
- Pump-to-PC Network Cable
- Power Supply