

Norm-Ject Syringes

Improved Design "Pharma" Grade – Sterile Packed Plastic Syringes

re-Usable



**NOT necessarily throw-a-way
NO "black tip" to clog or dry-out
NO Built-in obsolescence**

NORM-JECT® Syringes

NORM-JECT® syringes are latex-free, contain no rubber, no silicone oil, styrene or DEHP and are DNA-free. These syringes are the choice for any situation needing an inert, nonreactive syringe. Because of their composition, they are indicated for chromatography and many other laboratory procedures. They are more chemically resistant than rubber-tipped syringes and are manufactured from only laboratory grade polypropylene and polyethylene. These unique plastic syringes have a positive safety stop to prevent accidental spills. NORM-JECT® syringes are individually sterile strip packed.



Description	qty.	cat. #	18031
1ml Tuberculin ¹	100-pk.	22766	
3ml Luer Slip Centric Tip	100-pk.	22767	
5ml Luer Slip Centric Tip ²	100-pk.	22768	
10ml Luer Slip Eccentric Tip ²	100-pk.	22769	
20ml Luer Slip Eccentric Tip ²	100-pk.	22770	
30ml Luer Slip Eccentric Tip	50-pk.	22771	
50ml Luer Slip Eccentric Tip ²	30-pk.	22772	
3ml Luer Lock Tip	100-pk.	22773	
5ml Luer Lock Tip ¹	100-pk.	22774	
10ml Luer Lock Tip ²	100-pk.	22775	
20ml Luer Lock Tip ²	100-pk.	22776	
30ml Luer Lock Tip	50-pk.	22777	
50ml Luer Lock Tip ²	30-pk.	22778	



Great for use with Restek syringe filters! See page 363.

¹Dose saver design with low dead space plug on the piston to minimize waste.

²The 5ml has graduations to 6ml, 10ml has graduations to 12ml, 20ml has graduations to 24ml and 50ml has graduations to 60ml.