

## NUTECH 3602 SERIES AUTOSAMPLER

This instrument works with either the Nutech 8900 or Nutech 3551, it holds up to 16 canisters (all sizes) and Tedlar™ bags. The autosampler is a vertical design which saves valuable laboratory space and makes it easier to add and remove canisters. The canister hold plate is made in stainless steel and can be flipped up when not in use.



### FEATURES

- Space saving
- Easy to add or remove canisters
- Automated analysis of 16 Silonite® or SUMMA passivated canisters in accordance with EPA methods TO-14 and TO-15
- Back-flushing of each line after analysis under software control
- Automatic leak-check prior to opening canister valves
- Interfaces with 8900/3551 software

### SPECIFICATIONS

- Works with 8900 or 3551 Preconcentrator
- Two models available:
  - Model 3602 tower design
  - Model 3603 tower design (smaller footprint)
- 16 channels
- Leak check each position with QA/QC reports (controlled by preconcentrator software)
- Flush each position individually (controlled by preconcentrator software)
- 16 position valve heated with temperature control @  $\pm 2^{\circ}\text{C}$
- Dimensions:
  - Model 3602: 35" (W)  $\times$  59" (H)  $\times$  18" (D)
  - Model 3603: 30" (W)  $\times$  59" (H)  $\times$  15" (D)
- Weight:
  - Model 3602: 33 lbs.
  - Model 3603: 29 lbs.

